

LB102895

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

5/21/2019 13:02

Test: ~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	95.893	0.0	0.101	
ICB1	0.390	0.0	0.004	
CCV1	239.268	0.0	0.247	
CCB1	0.261	0.0	0.004	
PB119892BL	0.099	0.0	0.004	
K2641-12	0.074	0.0	0.004	
K2641-12DUP	0.255	0.0	0.004	
K2647-10	0.274	0.0	0.004	
K2913-04	0.184	0.0	0.004	
K2922-04	0.384	0.0	0.004	
K2922-08	-0.045	0.0	0.004	
PB119893BL	0.125	0.0	0.004	
K2913-06	0.346	0.0	0.004	
K2913-06DUP	0.580	0.0	0.005	
CCV2	239.641	0.0	0.247	
CCB2	0.459	0.0	0.004	
K2935-01	0.239	0.0	0.004	
PB119936BL	0.186	0.0	0.004	
K2951-16	0.191	0.0	0.004	
K2951-16DUP	0.223	0.0	0.004	
K2956-02	0.371	0.0	0.004	
PB119935BL	0.255	0.0	0.004	
K2951-04	0.261	0.0	0.004	
K2951-04DUP	0.176	0.0	0.004	
K2951-10	0.190	0.0	0.004	
K2951-14	0.267	0.0	0.004	
CCV3	242.048	0.0	0.250	
CCB3	0.383	0.0	0.004	

N	28
Mean	29.392
SD	76.5634
CV%	260.49

As

Aquakem v. 7.2AQ1

Results from time period:

Tue May 21 12:33:56 2019

Tue May 21 12:48:44 2019

Sample Id	Sam/Ctr/c/	Test short name	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	-0.2491	µg/l	5/21/2019 10:21:48	
5.OPPBCN	A	Reactive CN	P	4.5569	µg/l	5/21/2019 10:21:49	
10PPBCN	A	Reactive CN	P	9.3769	µg/l	5/21/2019 10:21:50	
50PPBCN	A	Reactive CN	P	55.8326	µg/l	5/21/2019 10:21:51	
100PPBCN	A	Reactive CN	P	96.5146	µg/l	5/21/2019 10:21:52	
250PPBCN	A	Reactive CN	P	247.6748	µg/l	5/21/2019 10:21:53	
500PPBCN	A	Reactive CN	P	501.2933	µg/l	5/21/2019 10:21:54	
ICV1	S	Reactive CN	P	95.8934	µg/l	5/21/2019 12:33:56	
ICB1	S	Reactive CN	P	0.3903	µg/l	5/21/2019 12:33:57	
CCV1	S	Reactive CN	P	239.2683	µg/l	5/21/2019 12:33:58	
CCB1	S	Reactive CN	P	0.2608	µg/l	5/21/2019 12:33:59	
PB119892BL	S	Reactive CN	P	0.099	µg/l	5/21/2019 12:34:00	
K2641-12	S	Reactive CN	P	0.0742	µg/l	5/21/2019 12:34:01	
K2641-12DUP	S	Reactive CN	P	0.255	µg/l	5/21/2019 12:34:02	
K2647-10	S	Reactive CN	P	0.2745	µg/l	5/21/2019 12:34:03	
K2913-04	S	Reactive CN	P	0.1842	µg/l	5/21/2019 12:34:04	
K2922-04	S	Reactive CN	P	0.3845	µg/l	5/21/2019 12:34:05	
K2922-08	S	Reactive CN	P	-0.0452	µg/l	5/21/2019 12:34:06	
PB119893BL	S	Reactive CN	P	0.1248	µg/l	5/21/2019 12:41:28	
K2913-06	S	Reactive CN	P	0.3464	µg/l	5/21/2019 12:41:29	
K2913-06DUP	S	Reactive CN	P	0.58	µg/l	5/21/2019 12:41:30	
CCV2	S	Reactive CN	P	239.6414	µg/l	5/21/2019 12:41:31	
CCB2	S	Reactive CN	P	0.4592	µg/l	5/21/2019 12:41:32	
K2935-01	S	Reactive CN	P	0.2393	µg/l	5/21/2019 12:41:33	
PB119936BL	S	Reactive CN	P	0.186	µg/l	5/21/2019 12:41:34	
K2951-16	S	Reactive CN	P	0.191	µg/l	5/21/2019 12:41:35	
K2951-16DUP	S	Reactive CN	P	0.2229	µg/l	5/21/2019 12:41:36	
K2956-02	S	Reactive CN	P	0.3712	µg/l	5/21/2019 12:41:37	
PB119935BL	S	Reactive CN	P	0.2551	µg/l	5/21/2019 12:41:38	
K2951-04	S	Reactive CN	P	0.2608	µg/l	5/21/2019 12:48:35	
K2951-04DUP	S	Reactive CN	P	0.1759	µg/l	5/21/2019 12:48:36	
K2951-10	S	Reactive CN	P	0.1896	µg/l	5/21/2019 12:48:37	
K2951-14	S	Reactive CN	P	0.2673	µg/l	5/21/2019 12:48:38	
CCV3	S	Reactive CN	P	242.0479	µg/l	5/21/2019 12:48:39	
CCB3	S	Reactive CN	P	0.3827	µg/l	5/21/2019 12:48:40	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

5/21/2019 10:27


Test Total CN

Accepted 5/21/2019 10:27

Factor 984.7
Bias 0.004

Coeff. of det. 0.999738

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.004	-0.2491	0.0000	
2	5.0PPBCN	0.009	4.5569	5.0000	
3	10PPBCN	0.013	9.3769	10.0000	
4	50PPBCN	0.061	55.8326	50.0000	
5	100PPBCN	0.102	96.5146	100.0000	
6	250PPBCN	0.255	247.6748	250.0000	
7	500PPBCN	0.513	501.2933	500.0000	