



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102900

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2642-05	OK-01-052119	61	1.15	9.7	8.91	90.8
K2753-01	CONCRETE-1	62	1.00	2.00	2.00	100.0
K2753-02	CONCRETE-2	63	1.00	2.00	2.00	100.0
K2753-03	CONCRETE-3	64	1.00	2.00	2.00	100.0
K2753-04	CONCRETE-4	65	1.00	2.00	2.00	100.0
K2753-05	CONCRETE-5	66	1.00	2.00	2.00	100.0
K2952-10	BP-VPB202-GW-643-645	1	1.16	9.98	2.8	18.6
K2952-11	BP-VPB202-GW-763-765	2	1.14	9.97	2.64	17.0
K2970-01	TP-11-S1	3	1.18	9.93	8.82	87.3
K2970-02	TP-11-S2	4	1.13	9.62	6.89	67.8
K2970-03	TP-11-S	5	1.17	9.55	7.79	79.0
K2970-05	TP-11-D1	6	1.15	9.97	6.49	60.5
K2970-06	TP-11-D2	7	1.12	9.9	7.99	78.2
K2970-07	TP-11-D	8	1.15	9.94	5.24	46.5
K2970-09	TP-10-S1	9	1.16	9.81	8.4	83.7
K2970-10	TP-10-S2	10	1.15	9.94	5.59	50.5
K2970-11	TP-10-D1	11	1.19	9.7	6.35	60.6
K2970-12	TP-10-D2	12	1.11	9.79	8.96	90.4
K2970-13	TP-9-S1	13	1.18	9.65	8.06	81.2
K2970-14	TP-9-S2	14	1.13	9.87	8.44	83.6
K2970-15	TP-9-S	15	1.16	9.77	8.57	86.1
K2971-01	SUP-2A-G-1-052019	16	1.12	6.63	6.04	89.3
K2971-02	SUP-2A-G-2-052019	17	1.15	5.1	3.04	47.8
K2971-03	SUP-2A-G-3-052019	18	1.12	6.4	2.94	34.5
K2971-04	SUP-2A-G-4-052019	19	1.14	5.63	3.35	49.2
K2971-05	SUP-2A-G-5-052019	20	1.13	7.6	3.69	39.6
K2971-06	SUP-2A-G-6-052019	21	1.16	6.49	4.81	68.5
K2971-07	SUP-2A-G-7-052019	22	1.11	6.11	4.6	69.8
K2971-08	SUP-2A-G-8-052019	23	1.16	6.05	4.14	60.9
K2971-09	SUP-2A-G-9-052019	24	1.11	6.03	5.66	92.5
K2971-10	SUP-2A-G-10-052019	25	1.19	6.5	5.91	88.9
K2971-11	SUP-2A-G-11-052019	26	1.18	6.54	6.08	91.4
K2971-12	SUP-2A-C-1-052019	27	1.14	9.75	8.87	89.8
K2971-13	SUP-2A-C-2-052019	28	1.19	9.58	8.86	91.4
K2971-14	SUP-2A-C-3-052019	29	1.15	9.84	9.04	90.8
K2971-15	SUP-2A-C-4-052019	30	1.1	9.58	8.9	92.0
K2971-16	SUP-2C-G-1-052019	31	1.18	5.84	5.47	92.1
K2971-17	SUP-2C-C-1-052019	32	1.15	9.52	8.84	91.9
K2972-01	DRUM-1-3	67	1.00	2.00	2.00	100.0



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102900

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2974-01	BP-VPB202-GW-823-825	33	1.19	9.94	2.87	19.2
K2976-01	EP1-SB-E2-3R(14-17)	34	1.13	4.45	3.9	83.4
K2976-02	EP1-SB-E2-3R(17-20)	35	1.13	5.82	5.14	85.5
K2976-03	EP1-SB-E3-4R(14-17)	36	1.11	5.39	4.79	86.0
K2976-04	EP1-SB-E3-4R(17-20)	37	1.18	5.17	4.5	83.2
K2976-05	EP1-SB-E-6R(14-17)	38	1.18	5.69	4.99	84.5
K2976-06	EP1-SB-E-6R(17-20)	39	1.18	6.51	5.93	89.1
K2976-07	EP1-SB-E-7R(14-17)	40	1.12	5.45	4.77	84.3
K2976-08	EP1-SB-E-7R(17-20)	41	1.15	5.39	4.76	85.1
K2976-09	EP1-SB-MW-ER(14-17)	42	1.15	9.9	8.77	87.1
K2976-10	K2976-09MS	43	1.15	9.9	8.77	87.1
K2976-11	K2976-09MSD	44	1.15	9.9	8.77	87.1
K2976-12	EP1-SB-MW-ER(17-20)	45	1.1	6.12	5.47	87.1
K2976-13	EP1-SB-MW-HR(14-17)	46	1.14	5.73	5.14	87.1
K2976-14	EP1-SB-MW-HR(17-20)	47	1.18	5.83	4.39	69.0
K2976-15	EP1-SB-MW-IR(14-17)	48	1.13	5.77	5.08	85.1
K2976-16	EP1-SB-MW-IR(17-20)	49	1.13	5.98	5.4	88.0
K2976-17	EP1-SB-SBR-2BR(14-17)	50	1.15	5.83	5.22	87.0
K2976-18	EP1-SB-SBR-2BR(17-20)	51	1.16	4.81	4.18	82.7
K2976-19	EP1-SB-SBR-3BR(14-17)	52	1.11	5.4	4.64	82.3
K2976-20	EP1-SB-SBR-3BR(17-20)	53	1.18	6.05	5.32	85.0
K2976-21	EP1-SBR-4BR(14-17)	54	1.17	6.09	5.2	81.9
K2976-22	EP1-SBR-4BR(17-20)	55	1.18	6.27	5.53	85.5
K2977-03	EP1-W-2R(14-17)	56	1.18	6.02	4.43	67.1
K2977-04	EP1-W-2R(17-20)	57	1.19	6.2	4.33	62.7
K2977-05	EP1-W-5R(14-17)	58	1.18	5.57	3.7	57.4
K2977-06	EP1-W-5R(17-20)	59	1.18	6.25	4.52	65.9
K2977-07	EP1-SB-FD-02	60	1.13	4.41	3.9	84.5
K2978-01	RO-SOIL	68	1.15	9.68	8.32	84.1
K2978-02	RO-SOIL	69	1.15	9.68	8.32	84.1
K2978-04	CONCRETE-COMP	70	1.00	2.00	2.00	100.0
K2978-06	TP-SOIL	71	1.11	9.73	7.17	70.3
K2978-07	TP-SOIL	72	1.11	9.73	7.17	70.3
K2983-01	SOIL-GRAVEL	73	1.00	2.00	2.00	100.0
K2983-03	OILY-DEBRIS	74	1.00	2.00	2.00	100.0
K2985-01	05-21-19-A	75	1.00	2.00	2.00	100.0
K2985-02	05-21-19-B	76	1.00	2.00	2.00	100.0
K2985-03	05-21-19-C	77	1.00	2.00	2.00	100.0
K2986-01	FESB-19(6.5-7)	78	1.13	9.65	8.14	82.3



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 05/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102900

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K2986-02	FESB-21 (4.5-5)	79	1.15	9.98	8.71	85.6
K2986-03	FEQC-052119	80	1.14	9.93	8.19	80.2
K2987-01	390516716	81	1.14	9.95	2.66	17.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$