



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/23/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 05/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:20
Out Date: 05/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102928

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2647-11	SU-03-052219-A	1	1.13	9.82	8.97	90.2
K2647-13	SU-03-052219-B	2	1.18	9.77	8.98	90.8
K2984-03	OILY-WATER	3	1.00	2.00	2.00	100.0
K2993-01	SB-1_2.0-2.5	4	1.13	9.79	8.18	81.4
K2993-02	SB-2_1.5-2.0	5	1.13	9.92	8.58	84.8
K2993-03	SB-3_1.5-2.0	6	1.17	9.81	8.65	86.6
K2993-04	SB-4_1.5-2.0	7	1.15	9.96	9.35	93.1
K2993-05	SB-5_1.5-2.0	8	1.16	9.62	8.67	88.8
K2993-06	SB-6_1.5-2.0	9	1.13	9.88	9.71	98.1
K2993-07	SB-7_1.0-1.5	10	1.14	9.75	9.31	94.9
K2993-08	SB-8_2.0-2.5	11	1.1	9.89	8.61	85.4
K2993-09	SB-9_1.0-1.5	12	1.18	9.93	8.76	86.6
K2993-10	SB-10_1.5-2.0	13	1.17	9.55	8.73	90.2
K2993-11	SB-11_2.0-2.5	14	1.16	9.96	8.83	87.2
K2993-12	SB-12_1.0-1.5	15	1.19	9.98	9.35	92.8
K2993-13	SB-13_2.0-2.5	16	1.18	9.96	9.02	89.3
K2993-14	SB-14_2.0-2.5	17	1.16	9.82	8.63	86.3
K2993-15	SB-15_1.5-2.0	18	1.13	9.69	8.52	86.3
K2993-16	COMP-1	19	1.12	9.62	8.43	86.0
K2993-17	COMP-2	20	1.11	9.72	8.63	87.3
K2993-18	COMP-3	21	1.13	9.81	8.76	87.9
K2994-01	TP-8-S1	22	1.13	9.6	6.58	64.3
K2994-02	TP-8-S2	23	1.14	9.64	7.48	74.6
K2994-03	TP-8-S	24	1.19	9.72	7.1	69.3
K2994-05	TP-8-D1	25	1.14	9.88	4.99	44.1
K2994-06	TP-8-D2	26	1.16	9.97	3.77	29.6
K2994-07	TP-8-D	27	1.14	9.79	5.63	51.9
K2994-09	TP-7-S1	28	1.18	9.97	6.23	57.5
K2994-10	TP-7-S2	29	1.13	9.55	5.84	55.9
K2994-11	TP-7-D1	30	1.16	9.96	5.18	45.7
K2994-12	TP-7-D2	31	1.16	9.57	5.59	52.7
K2994-13	TP-6-S1	32	1.19	9.8	6.86	65.9
K2994-14	TP-6-S2	33	1.13	9.83	6.68	63.8
K2994-15	TP-6-S	34	1.12	9.86	6.55	62.1
K2996-01	P-1 (9.5-10.0)	35	1.13	6.86	4.16	52.9
K2996-02	P-2 (0.5-1.0)	36	1.16	6.7	5.95	86.5
K2996-03	P-2 (4.5-5.0)	37	1.11	6.16	4.75	72.1
K2996-04	P-3 (0.5-1.0)	38	1.16	6.75	4.76	64.4
K2996-05	P-3 (4.5-5.0)	39	1.13	6.3	4.81	71.2



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/23/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 05/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:20
Out Date: 05/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102928

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2996-06	P-7 (0.5-1.0)	40	1.15	6.06	5.87	96.1
K2996-07	P-7 (4.5-5.0)	41	1.17	6.52	5.16	74.6
K2996-08	P-8 (9.5-10.0)	42	1.12	5.85	4.46	70.6
K2996-09	P-9 (0.5-1.0)	43	1.1	6.92	6.13	86.4
K2996-10	P-9 (4.5-5.0)	44	1.18	6.26	5.37	82.5
K2996-11	P-10 (0.5-1.0)	45	1.12	5.22	4.45	81.2
K2996-12	P-10 (4.5-5.0)	46	1.18	6.6	5.53	80.3
K2997-01	DPC-81626-FB-01	47	1.00	2.00	2.00	100.0
K2997-02	DPC-81626-FB-02	48	1.00	2.00	2.00	100.0
K2997-03	DPC-81626-FB-03	49	1.00	2.00	2.00	100.0
K2997-04	DPC-81626-FB-04	50	1.00	2.00	2.00	100.0
K2997-05	DPC-81626-FB-05	51	1.00	2.00	2.00	100.0
K2998-01	EWING-SP	52	1.16	9.83	8.27	82.0
K3002-01	COMP-Z1-20190418	53	1.19	9.71	8.27	83.1
K3002-02	COMP-Z2-20190419	54	1.15	9.96	7.5	72.1
K3002-03	COMP-Z3-20190420	55	1.14	9.69	7.41	73.3
K3002-04	21-C1-20190391	56	1.16	9.6	9.11	94.2
K3002-05	23-C3-20190399	57	1.13	9.67	7.25	71.7
K3002-06	22-C4-20190401	58	1.18	9.82	7.47	72.8
K3002-07	21-C5-20190409	59	1.14	9.64	5.93	56.4
K3002-08	22-C6-20190413	60	1.1	9.62	7.71	77.6
K3002-09	23-C7-20190417	61	1.19	9.77	6.09	57.1
K3004-01	CHRT-26835	62	1.17	9.82	7.78	76.4
K3008-01	COMP-19	63	1.13	9.57	6.47	63.3
K3008-02	VOC-19	64	1.14	5.25	4.58	83.7
K3008-03	COMP-20	65	1.16	9.96	8.56	84.1
K3008-04	VOC-20	66	1.18	6.59	5.77	84.8
K3008-05	COMP-21	67	1.13	9.83	7.61	74.5
K3008-06	VOC-21	68	1.11	6.66	6.05	89.0
K3010-01	FESB-20 (2.5-3)	69	1.15	9.93	7.78	75.5
K3010-02	FESB-23 (3-3.5)	70	1.19	9.54	8.02	81.8
K3010-03	FESB-22 (1.5-2)	71	1.14	9.92	7.93	77.3
K3010-04	k3010-02MS	72	1.19	9.54	8.02	81.8
K3010-05	k3010-02MSD	73	1.19	9.54	8.02	81.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$