



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/24/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102960

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2643-09	OR-03-052319-A	1	1.12	9.7	8.75	88.9
K2643-11	OR-03-052319-B	2	1.14	9.61	8.4	85.7
K3011-01	390506740	3	1.18	9.63	4.02	33.6
K3012-01	390508703	4	1.19	9.97	1.67	5.5
K3013-01	390508717.2	5	1.13	9.61	1.47	4.0
K3019-01	RT-3422	6	1.15	9.96	9.1	90.2
K3019-02	RT-3422	7	1.15	9.96	9.1	90.2
K3019-04	RBR200036	8	1.17	9.93	8.47	83.3
K3019-05	RBR200036	9	1.17	9.93	8.47	83.3
K3019-07	RT-1465	10	1.16	9.96	9.35	93.1
K3019-08	RT-1465	11	1.16	9.96	9.35	93.1
K3023-01	SSSI-0523-S-TCLP-11	12	1.16	9.77	7.66	75.5
K3026-01	TP-5-S1	13	1.11	9.88	8.36	82.7
K3026-02	TP-5-S2	14	1.18	9.58	8.27	84.4
K3026-03	TP-5-S	15	1.13	9.53	6.82	67.7
K3026-05	TP-5-D1	16	1.14	9.9	6.63	62.7
K3026-06	TP-5-D2	17	1.18	9.95	4.84	41.7
K3026-07	TP-5-D	18	1.14	9.71	6.76	65.6
K3026-09	TP-4-S1	19	1.16	9.7	6.92	67.4
K3026-10	TP-4-S2	20	1.18	9.82	5.77	53.1
K3026-11	TP-4-D1	21	1.1	9.53	6.33	62.0
K3026-12	TP-4-D2	22	1.14	9.98	6.85	64.6
K3026-13	TP-3-S1	23	1.14	9.95	7.84	76.0
K3026-14	TP-3-S2	24	1.15	9.94	7.59	73.3
K3026-15	TP-3-S	25	1.12	9.89	7.61	74.0
K3036-01	FESB-26 (2.5-3)	26	1.16	9.6	8.37	85.4
K3036-02	FEQC-DUP-052319	27	1.2	9.93	8.52	83.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$