



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 05/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103024

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3006-01	BFH86	1	1.13	9.71	4.63	40.8
K3006-02	BFH88	2	1.14	9.56	5.37	50.2
K3006-03	BFH88MS	3	1.14	9.56	5.37	50.2
K3006-04	BFH88MSD	4	1.14	9.56	5.37	50.2
K3006-05	BFHC0	5	1.16	9.68	5.17	47.1
K3006-06	BFHC1	6	1.12	9.63	4.49	39.6
K3006-07	BFHC4	7	1.15	9.95	6.13	56.6
K3006-08	BFHD0	8	1.18	9.96	4.7	40.1
K3006-09	VHBLK01	9	1.00	2.00	2.00	100.0
K3006-10	BFH78	10	1.13	9.98	6.83	64.4
K3006-11	BFH76	11	1.14	9.89	6.09	56.6
K3006-12	BFH74	12	1.18	9.95	5.02	43.8
K3006-13	BFH72	13	1.16	9.79	4.82	42.4
K3006-14	BFH80	14	1.19	9.87	7.03	67.3
K3006-15	BFH81	15	1.12	9.94	6.25	58.2
K3006-16	BFH82	16	1.18	9.89	7.28	70.0
K3006-17	BFH84	17	1.18	9.97	4.57	38.6
K3006-18	BFHB3	18	1.14	9.94	6.24	58.0
K3006-19	BFHC6	19	1.1	9.64	5.3	49.2
K3006-20	BFHC7	20	1.15	9.94	4.17	34.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$