



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/31/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:15
In Date: 05/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:37
Out Date: 05/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103089

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3100-01	ETMX2	1	1.14	9.62	7.45	74.4
K3100-02	ETMX3	2	1.15	9.59	8.28	84.5
K3100-03	ETMX4	3	1.14	9.7	9.15	93.6
K3100-04	ETMX4MS	4	1.14	9.7	9.15	93.6
K3100-05	ETMX4MSD	5	1.14	9.7	9.15	93.6
K3100-06	ETMX5	6	1.15	9.87	8.16	80.4
K3100-07	ETMX6	7	1.13	9.68	8.11	81.6
K3100-08	ETMX7	8	1.12	9.86	8.51	84.6
K3100-09	ETMX8	9	1.13	9.5	9.22	96.7
K3100-10	ETMX9	10	1.13	9.63	8.3	84.4
K3100-11	ETMY0	11	1.14	9.77	8.63	86.8
K3100-12	ETMY1	12	1.15	9.66	8.26	83.5
K3100-13	ETMY2	13	1.14	9.63	8.21	83.3
K3100-14	ETMY3	14	1.13	9.85	9.1	91.4
K3100-15	ETMY4	15	1.15	9.71	7.48	73.9
K3100-16	ETMY5	16	1.14	9.98	8.86	87.3
K3100-17	ETMY6	17	1.13	9.95	8.41	82.5
K3100-18	ETMY7	18	1.13	9.94	8.2	80.2
K3100-19	ETMY8	19	1.15	9.55	7.79	79.0
K3100-20	ETMZ0	20	1.14	9.78	7.84	77.5
K3100-21	ETMZ1	21	1.12	9.79	8.07	80.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$