

OVEN TEMP1: 103 °C TEMP: 105 °C

DATE & TIME IN: 06/04/2019 15:20 OUT: 06/05/2019 07:10

OVEN TEMP2: 105 °C TEMP: 105 °C

DATE & TIME IN: 06/05/2019 08:30 OUT: 06/05/2019 10:00

Run Number: LB103190

BalanceID: WC SC-3

OvenID: WC-OVEN-2

Dish #	Lab ID	Client ID	Empty Dish Weight (g)	Final Empty Dish Weight (g)	Sample Vol (ml)	1st Dish+Sample weight after Drying @103-@105°C (g)	2nd Dish+Sample weight after Drying @103-@105°C (g)	Final Dish+Sample weight after Drying @103-@105°C (g)	Weight (g)	Result (mg/L)
1	LB103190BL	LB103190BL	82.5564	82.5564	100	82.5564	82.5564	82.5564	0	0
2	K3182-07	MW-10	89.1618	89.1618	100	89.1916	89.1916	89.1916	0.0298	298
3	K3182-07DUP	MW-10DUP	75.5564	75.5564	100	75.586	75.586	75.586	0.0296	296

A = Sample Volume (ml)

B = Final Empty Dish Weight (g)

C = Final Dish+Sample weight after Drying @103-@105°C (g)

$$\text{Result mg/L} = ((C - B) / A) * 1000 * 1000$$

WORKLIST(Hardcopy Internal Chain)

WorkList Name : ts060419

WorkList ID : 126942

Date : 06-04-2019 14:22:07

XB/03190

Due Date	Matrix	Sample	Test	Preservative	Customer	Storage Location	Customer Sample	Collect Date	Method
06/11/2019	Water	K3182-07	TS	Cool 4 deg C	LIRO01	Q21	MW-10	06/04/2019	SM2540 B

Date/Time 06/04/19 14:40
 Received by: 
 Relinquished by:

Date/Time 06/04/19 15:30
 Received by: 
 Relinquished by: 