



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/6/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:25
In Date: 06/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:32
Out Date: 06/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103207

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3199-01	ETNB4	1	1.12	9.69	7.66	76.3
K3199-02	ETNB5	2	1.15	9.86	8.34	82.5
K3199-03	ETNB6	3	1.14	9.57	7.94	80.7
K3199-04	ETNB7	4	1.16	9.77	7.83	77.5
K3199-05	ETNB8	5	1.18	9.75	7.81	77.4
K3199-06	ETNB9	6	1.18	9.97	8.39	82.0
K3199-07	ETNC0	7	1.15	9.95	8.59	84.5
K3199-08	ETNC1	8	1.11	9.96	8.45	82.9
K3199-09	ETNC2	9	1.13	9.68	8.22	82.9
K3199-10	ETNC3	10	1.14	9.94	8.28	81.1
K3199-11	ETND4	11	1.17	9.71	8.31	83.6
K3199-12	ETND4MS	12	1.17	9.71	8.31	83.6
K3199-13	ETND4MSD	13	1.17	9.71	8.31	83.6
K3199-14	ETND5	14	1.18	9.87	8.04	78.9
K3199-15	ETND6	15	1.17	9.53	8.38	86.2
K3199-16	ETND7	16	1.1	9.92	8.74	86.6
K3199-17	ETND9	17	1.16	9.98	8.23	80.2
K3199-18	ETNE0	18	1.17	9.61	8.12	82.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$