



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/6/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:10
In Date: 06/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 06/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103209

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3201-01	C0AB1	1	1.14	9.64	9.06	93.2
K3201-02	C0AB3	2	1.19	9.78	9.31	94.5
K3201-03	C0AB4	3	1.12	9.69	9.3	95.4
K3201-04	C0AB5	4	1.14	9.82	9.23	93.2
K3201-05	C0AB6	5	1.13	9.9	9.29	93.0
K3201-06	C0AB7	6	1.15	9.63	8.38	85.3
K3201-07	C0AB8	7	1.11	9.55	7.94	80.9
K3201-08	C0AB9	8	1.13	9.58	7.68	77.5
K3201-09	C0AC0	9	1.18	9.55	8.43	86.6
K3201-10	C0AC1	10	1.11	9.73	7.59	75.2
K3201-11	C0AC2	11	1.17	9.62	8.42	85.8
K3201-12	C0AC3	12	1.16	9.83	8.75	87.5
K3201-13	C0AC4	13	1.12	9.68	9.02	92.3
K3201-14	C0AC5	14	1.11	9.65	7.35	73.1
K3201-15	C0AC6	15	1.15	9.71	7.62	75.6
K3201-16	C0AC7	16	1.13	9.51	7.54	76.5
K3201-17	C0AC8	17	1.18	9.71	8.23	82.6
K3201-18	C0AC9	18	1.15	9.55	8.39	86.2
K3201-19	C0AD0	19	1.18	9.6	8.5	86.9
K3201-20	C0AD1	20	1.16	9.72	8.41	84.7

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$