



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:35
In Date: 06/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:20
Out Date: 06/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103229

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3205-01	MH0AX1	1	1.15	9.66	7.4	73.4
K3205-02	MH0AX2	2	1.13	9.7	8.51	86.1
K3205-03	MH0AX3	3	1.16	9.73	8.74	88.4
K3205-04	MH0AX6	4	1.15	9.96	9.27	92.2
K3205-05	MH0AX8	5	1.17	9.59	8.07	81.9
K3205-06	MH0AX9	6	1.15	9.94	8.95	88.7
K3205-07	MH0AY0	7	1.12	9.53	8.73	90.5
K3205-08	MH0AY1	8	1.13	9.71	8.49	85.8
K3205-09	MH0AY2	9	1.17	9.51	8.72	90.5
K3205-10	MH0AY3	10	1.13	9.76	8.91	90.2
K3205-11	MH0AY4	11	1.17	9.91	9.18	91.6
K3205-12	MH0AY5	12	1.18	9.94	9.15	91.0
K3205-13	MH0AY6	13	1.18	9.83	8.66	86.5
K3205-14	MH0AY7	14	1.18	9.75	8.22	82.1
K3205-15	MH0AY8	15	1.1	9.6	8.11	82.5
K3205-16	MH0AY9	16	1.15	9.65	8.39	85.2
K3205-17	MH0AZ1	17	1.19	9.96	9.14	90.6
K3205-18	MH0AZ2	18	1.18	9.96	8.94	88.4
K3205-19	MH0AZ3	19	1.13	9.58	8.65	89.0
K3205-20	MH0AZ4	20	1.13	9.78	9.26	94.0
K3205-21	MH0AZ4D	21	1.13	9.78	9.26	94.0
K3205-22	MH0AZ4S	22	1.13	9.78	9.26	94.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$