



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:02
In Date: 06/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 06/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103270

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3231-01	DB8H2	1	1.00	2.00	2.00	100.0
K3231-02	DB8H3	2	1.00	2.00	2.00	100.0
K3231-03	DB8H4	3	1.00	2.00	2.00	100.0
K3231-04	DB8L3	4	1.15	9.92	8.94	88.8
K3231-05	DB8L4	5	1.14	9.56	8.43	86.6
K3231-06	DB8L4MS	6	1.14	9.56	8.43	86.6
K3231-07	DB8L4MSD	7	1.14	9.56	8.43	86.6
K3231-08	DB8L5	8	1.14	9.67	8.76	89.3
K3231-09	DB8L6	9	1.14	9.81	8.56	85.6
K3231-10	DB8L7	10	1.15	9.62	8.71	89.3
K3231-11	DB8L8	11	1.14	9.73	8.84	89.6
K3231-12	DB8L9	12	1.14	9.54	8.65	89.4
K3231-13	DB8M0	13	1.14	9.61	8.77	90.1
K3231-14	DB8M1	14	1.12	9.79	8.61	86.4
K3231-15	DB8M2	15	1.14	9.82	8.15	80.8
K3231-16	DB8M3	16	1.14	9.7	7.62	75.7
K3231-17	DB8M4	17	1.13	9.78	8.01	79.5
K3231-18	DB8M5	18	1.14	9.75	8.73	88.2
K3231-19	DB8M6	19	1.14	9.86	8.66	86.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$