



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/11/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:35
In Date: 06/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103310

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3252-04	TK11-CONTAINMENT	38	1.00	2.00	2.00	100.0
K3252-05	TK13-CONTAINMENT	39	1.00	2.00	2.00	100.0
K3252-06	TK7-CONTAINMENT	40	1.00	2.00	2.00	100.0
K3265-01	B-01-1.5-2.0	1	1.14	9.9	8.27	81.4
K3265-02	B-01-(0-37) COMP	2	1.16	9.97	8.53	83.7
K3265-03	B-04-2.5-3.0	3	1.15	9.74	8.43	84.7
K3265-04	B-04-COMP	4	1.14	9.94	8.62	85.0
K3265-05	B-03-2.5-3.0	5	1.14	9.98	7.05	66.9
K3265-06	B-03-COMP	6	1.12	9.74	7.66	75.9
K3265-07	B-02-2.5-3.0	7	1.16	9.68	8.43	85.3
K3265-08	B-02-COMP	8	1.12	9.91	8.58	84.9
K3271-01	BAT-B1	9	1.19	9.64	8.79	89.9
K3271-02	BAT-B1	10	1.19	9.94	9.03	89.6
K3271-03	BAT-B1-GEO	11	1.15	9.73	8.53	86.0
K3273-01	SB-6_1.0-1.5_060719	12	1.13	9.55	8.27	84.8
K3273-02	SB-6_12.5-13.0_060719	13	1.14	9.57	9.02	93.5
K3273-03	SB-1_1.0-1.5_060719	14	1.15	9.95	8.92	88.3
K3273-04	SB-1_9.0-9.5_060719	15	1.19	9.69	8.83	89.9
K3273-05	SB-4_1.0-1.5_060719	16	1.13	9.97	9.09	90.0
K3273-06	SB-4_9.5-10.0_060719	17	1.13	9.73	9.24	94.3
K3273-07	SB-5_1.0-1.5_060719	18	1.14	9.6	8.84	91.0
K3273-08	SB-5_9.5-10.0_060719	19	1.17	9.77	8.72	87.8
K3273-09	SB-2_1.0-1.5_060719	20	1.17	9.88	8.52	84.4
K3273-10	SB-2_9.5-10.0_060719	21	1.1	9.51	9.08	94.9
K3273-11	SB-3_1.0-1.5_060719	22	1.17	9.97	9.02	89.2
K3273-12	SB-3_9.0-9.5_060719	23	1.13	9.56	8.93	92.5
K3273-13	SB-8_1.0-1.5_060719	24	1.18	9.59	8.81	90.7
K3273-14	SB-8_9.5-10.0_060719	25	1.11	9.79	8.9	89.7
K3273-15	SB-7_1.0-1.5_060719	26	1.18	9.88	9.06	90.6
K3273-16	SB-7_9.5-10.0_060719	27	1.2	9.98	8.55	83.7
K3273-17	SB-9_0.5-1.0_060719	28	1.13	9.76	7.13	69.5
K3273-18	SB-9_1.0-1.5_060719	29	1.18	9.94	8.00	77.9
K3273-19	SB-10_0.5-1.0_060719	30	1.15	9.52	6.19	60.2
K3273-20	SB-10_1.0-1.5_060719	31	1.12	9.75	8.32	83.4
K3274-01	SB-11_0.5-1.0_060719	32	1.16	9.95	7.88	76.5
K3274-02	SB-11_1.0-1.5_060719	33	1.11	9.91	8.83	87.7
K3274-03	SB-12_0.5-1.0_060719	34	1.14	9.6	7.19	71.5
K3274-04	SB-12_1.0-1.5_060719	35	1.18	9.63	7.89	79.4
K3274-05	SB-13_0.5-1.0_060719	36	1.15	9.71	6.52	62.7



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 6/11/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:35
In Date: 06/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 06/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103310

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K3274-06	SB-13_2.0-2.5_060719	37	1.19	9.53	8.08	82.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$