



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: JIGNESH

Date: 6/21/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 17:25

In Date: 06/20/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:17

Out Date: 06/21/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB103535

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3158-05	OK-03-062019	38	1.12	9.71	8.94	91.0
K3445-01	EXT-008-SW116-061719	1	1.15	9.82	9.06	91.2
K3445-02	EXT-001-SW079-061719	2	1.15	9.54	8.89	92.3
K3445-03	EXT-008-SW117-061719	3	1.13	9.62	8.5	86.8
K3445-04	EXT-008-SW112-061719	4	1.11	9.72	9.09	92.7
K3445-05	EXT-008-SW113-061719	5	1.13	9.73	8.32	83.6
K3445-06	PST-007-SW050-061819	6	1.12	9.74	8.96	91.0
K3445-07	PST-007-B257-061819	7	1.14	9.61	8.77	90.1
K3445-08	PST-007-SW048-061819	8	1.13	9.86	9.22	92.7
K3445-09	PST-007-SW049-061819	9	1.14	9.53	8.98	93.4
K3445-10	EXT-007-SW124-061819	10	1.14	9.78	8.88	89.6
K3445-11	PST-029-B217-061719	11	1.12	9.63	8.32	84.6
K3445-12	PST-029-B052-061719	12	1.13	9.61	8.78	90.2
K3445-13	EXT-029-SW170-061719	13	1.12	9.98	9.32	92.6
K3445-14	EXT-029-SW171-061719	14	1.13	9.73	8.89	90.2
K3445-15	PST-001-B198-061719	15	1.12	9.58	8.99	93.0
K3445-16	PST-001-SW100-061719	16	1.12	9.64	8.7	89.0
K3445-17	EXT-001-SW077-061719	17	1.12	9.54	8.65	89.4
K3445-18	EXT-001-SW078-061719	18	1.15	9.66	9.13	93.8
K3445-19	EXT-008-SW122-061719	19	1.12	9.91	9.16	91.5
K3445-20	EXT-008-SW110-061719	20	1.14	9.8	8.93	90.0
K3445-20DUP	EXT-008-SW110-061719DUP	21	1.14	9.83	8.97	90.1
K3446-01	EXT-007-SW126-061819	22	1.13	9.93	8.93	88.6
K3446-02	PST-007-B255-061819	23	1.12	9.81	8.73	87.6
K3446-03	PST-007-B256-061819	24	1.12	9.84	9.05	90.9
K3446-04	PST-008-SW121-061819	25	1.12	9.61	9.18	94.9
K3446-05	PST-008-B251-061819	26	1.12	9.99	9.38	93.1
K3446-06	PST-008-B253-061819	27	1.12	9.98	9.24	91.6
K3446-07	PST-008-B254-061819	28	1.13	9.54	8.42	86.7
K3446-08	PST-008-B252-061819	29	1.14	9.99	9.16	90.6
K3446-09	PST-008-B036-061819	30	1.14	9.88	8.91	88.9
K3446-10	PST-008-SW122-061819	31	1.13	9.98	9.55	95.1
K3446-11	PST-008-SW124-061819	32	1.13	9.87	9.00	90.0
K3446-12	PST-008-SW123-061819	33	1.14	9.51	8.9	92.7
K3446-13	EXT-008-SW123-061819	34	1.12	9.93	8.79	87.1
K3446-14	EXT-008-SW123-061819-1	35	1.13	9.86	8.61	85.7
K3446-15	EXT-008-SW114-061819	36	1.13	9.61	9.02	93.0
K3446-15DUP	EXT-008-SW114-061819DUP	37	1.13	9.66	9.06	93.0
K3468-01	T-1	39	1.14	9.58	9.17	95.1



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: JIGNESH

Date: 6/21/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 17:25

In Date: 06/20/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104

Time OUT: 08:17

Out Date: 06/21/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB103535

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3468-02	T-2	40	1.13	9.71	9.22	94.3
K3468-03	T-3	41	1.18	9.68	9.25	94.9
K3468-04	T-4	42	1.1	9.71	9.26	94.8
K3468-05	T-5	43	1.16	9.6	9.06	93.6
K3468-06	T-6	44	1.1	9.8	9.04	91.3
K3470-01	SAMPLE-4	45	1.13	9.68	8.71	88.7
K3470-02	SAMPLE-4	46	1.13	9.68	8.71	88.7
K3470-03	SAMPLE-4A	47	1.13	9.68	8.71	88.7
K3470-05	SAMPLE-5	48	1.15	9.96	9.12	90.5
K3470-06	SAMPLE-5	49	1.15	9.96	9.12	90.5
K3470-07	SAMPLE-5A	50	1.15	9.96	9.12	90.5
K3470-09	SAMPLE-6	51	1.17	9.83	9.03	90.8
K3470-10	SAMPLE-6	52	1.17	9.83	9.03	90.8
K3470-11	SAMPLE-6A	53	1.17	9.83	9.03	90.8
K3471-01	WASTE	54	1.14	9.95	7.78	75.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$