



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/24/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 06/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:17
Out Date: 06/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103558

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2459-01	FK-BPM-04-002	1	1.15	9.79	8.88	89.5
K2459-02	FK-BPM-04-002	2	1.17	9.61	8.9	91.6
K3152-01	FK-BPM-04-003	3	1.15	9.55	8.63	89.0
K3152-02	FK-BPM-04-003	4	1.15	9.86	8.96	89.7
K3161-01	MATAUN-1	5	1.00	2.00	2.00	100.0
K3161-02	MATAUN-2	6	1.00	2.00	2.00	100.0
K3161-03	MATAUN-3	7	1.00	2.00	2.00	100.0
K3161-04	MATAUN-4	8	1.00	2.00	2.00	100.0
K3161-05	MATAUN-5	9	1.00	2.00	2.00	100.0
K3161-06	MATAUN-6	10	1.00	2.00	2.00	100.0
K3161-07	MATAUN-7	11	1.00	2.00	2.00	100.0
K3161-08	MATAUN-8	12	1.00	2.00	2.00	100.0
K3161-09	MATAUN-9	13	1.00	2.00	2.00	100.0
K3161-10	MATAUN-10	14	1.00	2.00	2.00	100.0
K3161-11	MATAUN-11	15	1.00	2.00	2.00	100.0
K3161-12	MATAUN-12	16	1.00	2.00	2.00	100.0
K3161-13	MATAUN-13	17	1.00	2.00	2.00	100.0
K3161-14	MATAUN-14	18	1.00	2.00	2.00	100.0
K3161-15	MATAUN-15	19	1.00	2.00	2.00	100.0
K3161-16	MATAUN-16	20	1.00	2.00	2.00	100.0
K3161-17	MATAUN-17	21	1.00	2.00	2.00	100.0
K3161-18	MATAUN-18	22	1.00	2.00	2.00	100.0
K3161-19	MATAUN-19	23	1.00	2.00	2.00	100.0
K3161-20	MATAUN-20	24	1.00	2.00	2.00	100.0
K3162-01	MATAUN-21	25	1.00	2.00	2.00	100.0
K3162-02	MATAUN-22	26	1.00	2.00	2.00	100.0
K3162-03	MATAUN-23	27	1.00	2.00	2.00	100.0
K3162-04	MATAUN-24	28	1.00	2.00	2.00	100.0
K3162-05	MATAUN-25	29	1.00	2.00	2.00	100.0
K3162-06	MATAUN-26	30	1.00	2.00	2.00	100.0
K3162-07	MATAUN-27	31	1.00	2.00	2.00	100.0
K3162-08	MATAUN-28	32	1.00	2.00	2.00	100.0
K3162-09	MATAUN-29	33	1.00	2.00	2.00	100.0
K3162-10	MATAUN-30	34	1.00	2.00	2.00	100.0
K3162-11	MATAUN-31	35	1.00	2.00	2.00	100.0
K3162-12	MATAUN-32	36	1.00	2.00	2.00	100.0
K3162-13	MATAUN-33	37	1.00	2.00	2.00	100.0
K3162-14	MATAUN-34	38	1.00	2.00	2.00	100.0
K3162-15	MATAUN-35	39	1.00	2.00	2.00	100.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/24/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 06/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:17
Out Date: 06/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103558

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3162-16	MATAUN-36	40	1.00	2.00	2.00	100.0
K3162-17	MATAUN-37	41	1.00	2.00	2.00	100.0
K3162-18	MATAUN-38	42	1.00	2.00	2.00	100.0
K3352-02	20190612-LOC-11	81	1.17	9.92	6.71	63.3
K3352-04	20190612-LOC-12	82	1.16	9.97	5.38	47.9
K3472-01	TP-5-1	43	1.14	9.62	8.07	81.7
K3472-02	TP-8	44	1.13	9.88	8.12	79.9
K3472-03	TP-8-1	45	1.14	9.6	7.97	80.7
K3472-04	TP-8-2	46	1.14	9.65	7.89	79.3
K3472-05	TP-8-3	47	1.14	9.95	7.81	75.7
K3475-01	SOILS-1	48	1.14	9.98	8.82	86.9
K3475-02	SOILS-2	49	1.14	9.83	8.33	82.7
K3475-03	SOILS-3	50	1.14	9.52	8.33	85.8
K3475-04	SOILS-4	51	1.13	9.93	8.48	83.5
K3477-05	SVOC-GPC-BLANK	52	1.00	2.00	2.00	100.0
K3477-06	PEST-GPC-BLANK	53	1.00	2.00	2.00	100.0
K3477-07	PEST-GPC-BLANK-SPIKE	54	1.00	2.00	2.00	100.0
K3477-08	PCB-GPC-BLANK	55	1.00	2.00	2.00	100.0
K3477-09	PCB-GPC-BLANK-SPIKE	56	1.00	2.00	2.00	100.0
K3477-10	SVOC-GPC2-BLANK	57	1.00	2.00	2.00	100.0
K3477-11	PEST-GPC2-BLANK	58	1.00	2.00	2.00	100.0
K3477-12	PEST-GPC2-BLANK-SPIKE	59	1.00	2.00	2.00	100.0
K3477-13	PCB-GPC2-BLANK	60	1.00	2.00	2.00	100.0
K3477-14	PCB-GPC2-BLANK -SPIKE	61	1.00	2.00	2.00	100.0
K3480-01	TCMS-WIPE-001-01	62	1.00	2.00	2.00	100.0
K3480-02	TCMS-WIPE-002-01	63	1.00	2.00	2.00	100.0
K3480-03	TCMS-WIPE-003-01	64	1.00	2.00	2.00	100.0
K3480-04	TCMS-WIPE-004-01	65	1.00	2.00	2.00	100.0
K3480-05	TCMS-WIPE-005-01	66	1.00	2.00	2.00	100.0
K3480-06	TCMS-WIPE-006-03	67	1.00	2.00	2.00	100.0
K3482-01	TP-1	68	1.13	9.82	8.73	87.5
K3482-02	TP-1	69	1.13	9.82	8.73	87.5
K3482-03	TP-2	70	1.13	9.69	8.88	90.5
K3482-04	TP-2	71	1.12	9.64	8.88	91.1
K3482-05	TP-3	72	1.14	9.51	8.26	85.1
K3482-06	TP-3	73	1.14	9.51	8.26	85.1
K3482-07	TP-4	74	1.14	9.85	8.82	88.2
K3482-08	TP-4	75	1.14	9.85	8.82	88.2
K3487-01	FK-PS-01-048	76	1.13	9.54	7.6	76.9



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/24/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 06/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:17
Out Date: 06/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103558

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3487-02	FK-PS-01-048	77	1.13	9.69	7.6	75.6
K3490-01	FK-PS-01-049	78	1.13	9.8	8.02	79.5
K3490-02	FK-PS-01-049	79	1.15	9.98	8.24	80.3
K3493-01	390620142	80	1.12	9.9	2.96	21.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$