



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 06/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103665

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3158-07	OK-02-062719	115	1.18	9.86	9.34	94.0
K3163-13	SU-01-062719-A	116	1.17	9.59	8.85	91.2
K3163-15	SU-01-062719-B	117	1.16	9.95	9.16	91.0
K3553-01	PE-4A	1	1.15	9.96	8.27	80.8
K3553-02	PE-4B	2	1.15	9.74	8.72	88.1
K3553-03	PE-5A	3	1.16	9.63	7.55	75.4
K3553-04	PE-5B	4	1.15	9.95	9.17	91.1
K3553-05	PE-6	5	1.13	9.55	7.45	75.1
K3553-06	PE-7A	6	1.14	9.82	8.38	83.4
K3553-07	PE-7B	7	1.13	9.82	7.86	77.4
K3553-08	PE-8	8	1.14	9.93	8.00	78.0
K3553-09	PE-8-1	9	1.15	9.98	8.71	85.6
K3553-10	PE-9A	10	1.14	9.79	8.99	90.8
K3553-11	PE-9B	11	1.13	9.53	7.52	76.1
K3553-12	PE-10A	12	1.12	9.55	7.94	80.9
K3553-13	PE-10B	13	1.13	9.81	7.65	75.1
K3553-14	DUPLICATE	14	1.12	9.72	8.43	85.0
K3554-01	SAMPLE-11	63	1.12	9.64	8.09	81.8
K3554-02	SAMPLE-11	64	1.12	9.64	8.09	81.8
K3554-03	SAMPLE-11A	65	1.12	9.64	8.09	81.8
K3554-05	SAMPLE-15	66	1.12	9.85	8.74	87.3
K3554-06	SAMPLE-15	67	1.12	9.85	8.74	87.3
K3554-07	SAMPLE-15A	68	1.12	9.85	8.74	87.3
K3554-09	SAMPLE-16	69	1.17	9.66	8.7	88.7
K3554-10	SAMPLE-16	70	1.17	9.66	8.7	88.7
K3554-11	SAMPLE-16A	71	1.17	9.66	8.7	88.7
K3555-01	SAMPLE-17	72	1.18	9.96	8.75	86.2
K3555-02	SAMPLE-17	73	1.18	9.96	8.75	86.2
K3555-03	SAMPLE-17A	74	1.18	9.96	8.75	86.2
K3555-05	SAMPLE-18	75	1.15	9.84	8.67	86.5
K3555-06	SAMPLE-18	76	1.15	9.84	8.67	86.5
K3555-07	SAMPLE-18A	77	1.15	9.84	8.67	86.5
K3555-09	SAMPLE-9	78	1.13	9.59	8.52	87.4
K3555-10	SAMPLE-9	79	1.13	9.59	8.52	87.4
K3555-11	SAMPLE-9A	80	1.13	9.59	8.52	87.4
K3555-13	SAMPLE-12	81	1.11	9.86	8.66	86.3
K3555-14	SAMPLE-12	82	1.11	9.86	8.66	86.3
K3555-15	SAMPLE-12A	83	1.11	9.86	8.66	86.3
K3555-17	SAMPLE-13	84	1.15	9.7	8.42	85.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 06/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103665

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3555-18	SAMPLE-13	85	1.15	9.7	8.42	85.0
K3555-19	SAMPLE-13A	86	1.15	9.7	8.42	85.0
K3555-21	SAMPLE-14	87	1.14	9.98	8.91	87.9
K3555-22	SAMPLE-14	88	1.14	9.98	8.91	87.9
K3555-23	SAMPLE-14A	89	1.14	9.98	8.91	87.9
K3558-01	EXT-015-SW084-062519	15	1.13	9.98	8.77	86.3
K3558-02	EXT-015-SW086-062519	16	1.12	9.64	8.95	91.9
K3558-03	EXT-015-SW087-062519	17	1.13	9.91	9.09	90.7
K3558-04	PST-015-B229-062519	18	1.11	9.96	9.06	89.8
K3558-05	EXT-015-SW085-062519	19	1.13	9.97	9.41	93.7
K3558-06	EXT-015-SW082-062519	20	1.14	9.95	9.00	89.2
K3558-07	EXT-015-SW081-062519	21	1.13	9.79	8.73	87.8
K3558-08	PST-032-B204-062519	22	1.13	9.96	9.36	93.2
K3558-09	EXT-032-SW173-062519	23	1.13	9.83	8.82	88.4
K3558-10	EXT-032-SW172-062519	24	1.13	9.85	8.86	88.6
K3558-11	EXT-032-SW174-062519	25	1.12	9.98	9.32	92.6
K3558-12	PST-032-B205-062519	26	1.13	9.79	9.11	92.1
K3558-13	PST-016-B237-062619	27	1.13	9.53	8.71	90.2
K3558-14	PST-016-B236-062619	28	1.12	9.97	9.47	94.4
K3558-15	PST-016-B235-062619	29	1.13	9.93	9.21	91.8
K3558-16	PST-016-B056-062619	30	1.12	9.92	9.42	94.3
K3558-17	EXT-016-SW100-062619	31	1.13	9.8	8.74	87.8
K3558-18	EXT-016-SW096-062619	32	1.13	9.65	8.69	88.7
K3558-19	EXT-016-SW095-062619	33	1.12	9.91	8.56	84.6
K3560-01	EXT-009-SW149-062419	34	1.13	9.95	9.26	92.2
K3560-03	EXT-009-SW147-062419	35	1.12	9.67	8.8	89.8
K3560-04	EXT-009-SW145-062419	36	1.13	9.67	8.68	88.4
K3560-05	EXT-009-SW142-062419	37	1.12	9.86	9.37	94.4
K3560-06	EXT-009-SW143-062419	38	1.14	9.92	9.31	93.1
K3560-07	EXT-009-SW137-062419	39	1.12	9.67	9.06	92.9
K3560-08	EXT-009-SW146-062419	40	1.13	9.7	8.88	90.4
K3560-09	PST-009-B221-062419	41	1.13	9.63	9.11	93.9
K3560-10	PST-009-B220-062419	42	1.13	9.61	9.03	93.2
K3560-11	PST-009-B219-062419	43	1.12	9.6	9.05	93.5
K3560-12	PST-009-B218-062419	44	1.13	9.65	9.07	93.2
K3560-13	EXT-009-SW134-062419	45	1.14	9.97	9.26	92.0
K3560-14	EXT-009-SW133-062419	46	1.15	9.85	8.98	90.0
K3560-15	PST-009-B223-062419	47	1.14	9.97	9.51	94.8
K3560-16	PST-009-B222-062419	48	1.14	9.97	9.16	90.8



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 06/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103665

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3560-17	EXT-012-SW153-062419	49	1.1	9.6	8.7	89.4
K3560-18	EXT-012-SW139-062419	50	1.15	9.91	8.9	88.5
K3560-19	EXT-011-SW138-062419	51	1.13	9.72	8.83	89.6
K3560-20	PST-011-B224-062419	52	1.13	9.59	8.05	81.8
K3561-01	PST-012-B226-062419	53	1.13	9.7	8.6	87.2
K3561-02	PST-012-B227-062419	54	1.13	9.81	7.7	75.7
K3561-03	PST-016-SW117-062519	55	1.14	9.79	8.9	89.7
K3561-04	PST-016-SW117-062519-1	56	1.13	9.68	8.58	87.1
K3561-05	EXT-016-SW099-062519	57	1.13	9.62	8.84	90.8
K3561-06	EXT-016-SW098-062519	58	1.15	9.93	8.65	85.4
K3561-07	PST-015-B228-062519	59	1.13	9.8	8.96	90.3
K3561-08	PST-015-SW109-062519	60	1.14	9.67	7.97	80.1
K3561-09	PST-015-B053-062519	61	1.13	9.6	8.52	87.2
K3561-10	EXT-016-SW097-062519	62	1.13	9.56	8.52	87.7
K3563-01	SOIL-COMP-1	90	1.17	9.96	9.03	89.4
K3563-02	SOIL-COMP-1	91	1.17	9.96	9.03	89.4
K3563-03	SOIL-COMP-2	92	1.11	9.91	8.33	82.0
K3563-04	SOIL-COMP-2	93	1.11	9.91	8.33	82.0
K3573-01	SR-SOIL-COMP-1	94	1.19	9.81	9.06	91.3
K3573-02	SR-SOIL-COMP-2	95	1.13	9.58	8.82	91.0
K3576-01	SP-1	96	1.12	9.6	9.23	95.6
K3576-02	SP-1	97	1.18	9.58	9.04	93.6
K3580-01	SB-11-9.5-10.0	118	1.14	9.93	8.04	78.5
K3580-02	SB-11-COMP	119	1.15	9.71	8.24	82.8
K3580-03	SB-12-9.5-10.0	120	1.19	9.58	8.36	85.5
K3580-04	SB-12-COMP	121	1.18	9.58	8.14	82.9
K3580-05	SB-13-9.5-10.0	122	1.12	9.59	8.25	84.2
K3580-06	SB-13-COMP	123	1.18	9.71	8.69	88.0
K3580-07	SB-14-9.5-10.0	124	1.15	9.76	9.03	91.5
K3580-08	SB-14-COMP	125	1.17	9.97	8.36	81.7
K3580-09	SB-15-9.5-10.0	126	1.13	9.89	9.18	91.9
K3580-10	SB-15-COMP	127	1.13	9.92	8.62	85.2
K3581-01	PE-11	98	1.18	9.96	7.9	76.5
K3581-02	PE-12B	99	1.1	9.51	7.78	79.4
K3581-03	PE-13B	100	1.13	9.65	7.74	77.6
K3581-04	PE-12A	101	1.14	9.91	8.22	80.7
K3581-05	PE-13A	102	1.15	9.52	8.37	86.3
K3581-06	PE-15B	103	1.14	9.5	7.99	81.9
K3581-07	PE-15A	104	1.18	9.92	9.31	93.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 06/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 06/28/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103665

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3581-08	PE-14	105	1.16	9.78	7.62	74.9
K3581-09	PE-17B	106	1.17	9.98	9.4	93.4
K3581-10	PE-16	107	1.14	9.76	8.46	84.9
K3581-11	PE-18B	108	1.15	9.82	8.14	80.6
K3581-12	PE-19B	109	1.11	9.92	8.73	86.5
K3581-13	PE-19A	110	1.11	9.78	8.6	86.4
K3581-14	PE-18A	111	1.16	9.74	8.86	89.7
K3581-15	DUPLICATE	112	1.17	9.52	8.3	85.4
K3581-16	PE-17A	113	1.13	9.57	8.28	84.7
K3582-01	FK-PS-01-049	114	1.13	9.8	8.02	79.5
K3583-01	01-A-01-B-01-C	128	1.00	2.00	2.00	100.0
K3583-02	02-A-02-B-02-C	129	1.00	2.00	2.00	100.0
K3584-01	FK-0D015-03-009	130	1.17	9.53	9.36	98.0
K3584-02	FK-0D015-03-009	131	1.18	9.66	9.52	98.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$