

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

6/27/2019 13:08

Test: *Reactive*
~~Total~~ CN

Sample Id	Result	Dil. 1 +	Response	Errors
ICV1	98.020	0.0	0.100	
ICB1	-0.517	0.0	0.002	
CCV1	246.048	0.0	0.247	
CCB1	-0.626	0.0	0.002	
PB120962BL	-0.789	0.0	0.002	
K3157-12	-0.807	0.0	0.002	
K3157-12DUP	-0.860	0.0	0.001	
K3157-14	-0.944	0.0	0.001	
K3157-16	-0.960	0.0	0.001	
K3157-18	-0.800	0.0	0.002	
K3157-20	-0.745	0.0	0.002	
K3159-10	-0.892	0.0	0.001	
K3159-12	-0.868	0.0	0.001	
K3507-02	-0.451	0.0	0.002	
CCV2	247.007	0.0	0.248	
CCB2	-0.600	0.0	0.002	
K3507-04	-0.879	0.0	0.001	
K3509-02	-0.868	0.0	0.001	
K3509-04	-0.921	0.0	0.001	
K3544-01	-0.840	0.0	0.001	
CCV3	246.052	0.0	0.247	
CCB3	-0.975	0.0	0.001	

N	22
Mean	37.399
SD	87.5436
CV%	234.08

Aquakem v. 7.2AQ1

Results from time period:

Thu Jun 27 11:20:15 2019

Thu Jun 27 13:08:39 2019

AS

Sample Id	Sam/Ctr/c/	Test short nan	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	-1.1631	µg/l	6/27/2019 11:44:59	
5.OPPBCN	A	Reactive CN	P	2.8814	µg/l	6/27/2019 11:45:00	
10PPBCN	A	Reactive CN	P	9.2152	µg/l	6/27/2019 11:45:01	
50PPBCN	A	Reactive CN	P	49.7543	µg/l	6/27/2019 11:45:02	
100PPBCN	A	Reactive CN	P	103.0523	µg/l	6/27/2019 11:45:03	
250PPBCN	A	Reactive CN	P	253.6177	µg/l	6/27/2019 11:45:04	
500PPBCN	A	Reactive CN	P	497.6422	µg/l	6/27/2019 11:45:05	
ICV1	S	Reactive CN	P	98.0198	µg/l	6/27/2019 12:58:16	
ICB1	S	Reactive CN	P	-0.5169	µg/l	6/27/2019 12:58:17	
CCV1	S	Reactive CN	P	246.0481	µg/l	6/27/2019 12:58:18	
CCB1	S	Reactive CN	P	-0.6257	µg/l	6/27/2019 12:58:19	
PB120962BL	S	Reactive CN	P	-0.7889	µg/l	6/27/2019 12:58:20	
K3157-12	S	Reactive CN	P	-0.807	µg/l	6/27/2019 12:58:21	
K3157-12DUP	S	Reactive CN	P	-0.8597	µg/l	6/27/2019 12:58:22	
K3157-14	S	Reactive CN	P	-0.9435	µg/l	6/27/2019 12:58:23	
K3157-16	S	Reactive CN	P	-0.9599	µg/l	6/27/2019 12:58:24	
K3157-18	S	Reactive CN	P	-0.7998	µg/l	6/27/2019 12:58:25	
K3157-20	S	Reactive CN	P	-0.745	µg/l	6/27/2019 12:58:26	
K3159-10	S	Reactive CN	P	-0.8924	µg/l	6/27/2019 13:05:51	
K3159-12	S	Reactive CN	P	-0.8681	µg/l	6/27/2019 13:05:52	
K3507-02	S	Reactive CN	P	-0.4505	µg/l	6/27/2019 13:05:53	
CCV2	S	Reactive CN	P	247.0069	µg/l	6/27/2019 13:05:54	
CCB2	S	Reactive CN	P	-0.5998	µg/l	6/27/2019 13:05:55	
K3507-04	S	Reactive CN	P	-0.8785	µg/l	6/27/2019 13:05:56	
K3509-02	S	Reactive CN	P	-0.8682	µg/l	6/27/2019 13:05:57	
K3509-04	S	Reactive CN	P	-0.9207	µg/l	6/27/2019 13:05:58	
K3544-01	S	Reactive CN	P	-0.8398	µg/l	6/27/2019 13:05:59	
CCV3	S	Reactive CN	P	246.0522	µg/l	6/27/2019 13:06:01	
CCB3	S	Reactive CN	P	-0.9755	µg/l	6/27/2019 13:08:39	

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

6/27/2019 11:53


Test Total CN

Accepted 6/27/2019 11:53

Factor 1005
Bias 0.002

Coeff. of det. 0.999832

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.001	-1.1631	0.0000	
2	5.0PPBCN	0.005	2.8814	5.0000	
3	10PPBCN	0.011	9.2152	10.0000	
4	50PPBCN	0.052	49.7543	50.0000	
5	100PPBCN	0.105	103.0523	100.0000	
6	250PPBCN	0.255	253.6177	250.0000	
7	500PPBCN	0.497	497.6422	500.0000	