



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/3/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:10
In Date: 07/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 07/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103777

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3550-05	PE-1A	1	1.13	9.82	8.51	84.9
K3550-06	PE-3A	2	1.16	9.72	8.66	87.6
K3618-01	HD-01-070219	3	1.14	9.82	9.26	93.5
K3637-01	RO-276	4	1.00	2.00	2.00	100.0
K3645-01	FK-0D015-03-010	5	1.15	9.65	9.52	98.5
K3645-02	FK-0D015-03-010	6	1.17	9.73	9.58	98.2
K3646-01	FK-0D015-03-011	7	1.12	9.55	9.39	98.1
K3646-02	FK-0D015-03-011	8	1.11	9.78	9.61	98.0
K3648-01	C-1-CORE-HOLE-1	9	1.16	9.88	8.14	80.0
K3648-02	C-2-CORE-HOLE-2	10	1.18	9.97	8.69	85.4
K3649-01	D-1-ROLL-OFF-SOCK-LILE	11	1.18	9.95	8.3	81.2
K3649-02	D-1X-ROLL-OFF-SOCK-LILE	12	1.18	9.95	8.3	81.2
K3650-01	FK-PS-01-051	13	1.1	9.63	7.74	77.8
K3650-02	FK-PS-01-051	14	1.18	9.95	7.81	75.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$