



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/8/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 07/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:10
Out Date: 07/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103807

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3508-01	OILY-SOIL-DRUM	1	1.14	9.62	9.15	94.5
K3622-01	OR-02-070319-A	2	1.15	9.66	9.22	94.8
K3622-03	OR-02-070319-B	3	1.18	9.81	9.36	94.8
K3624-01	PL-01-070319-A	4	1.15	9.96	9.08	90.0
K3624-03	PL-01-070319-B	5	1.14	9.58	8.74	90.0
K3624-05	PL-01-070319-C	6	1.17	9.62	8.72	89.3
K3624-07	PL-01-070319-D	7	1.12	9.68	8.76	89.3
K3653-01	HL6-2 (2.5-3)	8	1.18	9.81	8.97	90.3
K3653-02	HL6-1 (3.5-4)	9	1.15	9.5	8.6	89.2
K3653-03	HL5-2 (2.5-3)	10	1.16	9.65	8.89	91.0
K3653-04	HL5-1 (3-3.5)	11	1.14	9.55	8.49	87.4
K3653-05	HL4-2 (3-3.5)	12	1.15	9.89	8.86	88.2
K3653-06	HL4-1 (3-3.5)	13	1.16	9.97	8.91	88.0
K3653-07	HL3-2 (2.5-3)	14	1.13	9.82	8.77	87.9
K3653-08	HL3-1 (2.5-3)	15	1.18	9.66	8.39	85.0
K3653-09	HL2-2 (2.5-3)	16	1.1	9.93	9.24	92.2
K3653-10	HL2-1 (2.5-3)	17	1.15	9.97	9.08	89.9
K3653-11	HL1-2 (3-3.5)	18	1.12	9.63	8.37	85.2
K3653-12	HL1-1 (2.5-3)	19	1.16	9.6	8.65	88.7
K3654-01	390626049	20	1.18	9.98	1.38	2.3
K3661-01	390626716.1	21	1.17	9.96	1.47	3.4
K3663-01	VNJ-213	22	1.14	9.81	8.26	82.1
K3663-02	339	23	1.1	9.91	8.66	85.8
K3663-03	VNJ-213-339	24	1.13	9.55	7.98	81.4
K3668-01	TRANSFORMER-1	25	1.15	9.66	9.21	94.7
K3668-02	TRANSFORMER-1-A	26	1.15	9.66	9.21	94.7
K3668-03	TRANSFORMER-1-B	27	1.15	9.66	9.21	94.7
K3668-04	TRANSFORMER-1	28	1.15	9.66	9.21	94.7
K3668-05	SHORT-WALL-AREA	29	1.19	9.83	9.34	94.3
K3668-06	SHORT-WALL-AREA-A	30	1.19	9.83	9.34	94.3
K3668-07	SHORT-WALL-AREA-B	31	1.19	9.83	9.34	94.3
K3668-08	SHORT-WALL-AREA	32	1.19	9.83	9.34	94.3
K3668-09	TRANSFORMER-T2	33	1.17	9.98	9.51	94.7
K3668-10	TRANSFORMER-T2-A	34	1.17	9.98	9.51	94.7
K3668-11	TRANSFORMER-T2-B	35	1.17	9.98	9.51	94.7
K3668-12	TRANSFORMER-T2	36	1.17	9.98	9.51	94.7
K3677-01	OIL-DRUM-A	37	1.00	2.00	2.00	100.0
K3677-02	OIL-DRUM-B	38	1.00	2.00	2.00	100.0
K3677-03	OILY-DEBRIS	39	1.00	2.00	2.00	100.0



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/8/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 07/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:10
Out Date: 07/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB103807

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3681-01	PILEA-07-01-19	40	1.15	9.97	9.04	89.5
K3681-02	TOP-SOIL-07-01-19	41	1.15	9.95	9.38	93.5
K3681-03	STOCK-PILE-07-01-19	42	1.19	9.92	9.29	92.8
K3682-01	D-2-ROLL-OFF-2-STOCK-PI	43	1.16	9.93	7.4	71.2
K3682-02	D-2X-ROLL-OFF-2-STOCK-P	44	1.14	9.86	7.31	70.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$