



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:42
In Date: 07/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-Oven

QC:LB103941

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3795-01	MH11K8	1	1.00	2.00	2.00	100.0
K3795-02	MH11K8D	2	1.00	2.00	2.00	100.0
K3795-03	MH11K8S	3	1.00	2.00	2.00	100.0
K3795-04	MH11N4	4	1.00	2.00	2.00	100.0
K3795-05	MH11N5	5	1.00	2.00	2.00	100.0
K3795-06	MH11N6	6	1.00	2.00	2.00	100.0
K3795-07	MH11N7	7	1.00	2.00	2.00	100.0
K3795-08	MH11N8	8	1.00	2.00	2.00	100.0
K3795-09	MH11N9	9	1.00	2.00	2.00	100.0
K3795-10	MH11P0	10	1.00	2.00	2.00	100.0
K3795-11	MH11P1	11	1.00	2.00	2.00	100.0
K3795-12	MH11P2	12	1.00	2.00	2.00	100.0
K3795-13	MH11P3	13	1.00	2.00	2.00	100.0
K3795-14	MH11P4	14	1.00	2.00	2.00	100.0
K3795-15	MH11P5	15	1.00	2.00	2.00	100.0
K3795-16	MH11P6	16	1.00	2.00	2.00	100.0
K3795-17	MH11P7	17	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$