



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:25
In Date: 07/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 07/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB103945

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3620-01	CL-01-071219-A	1	1.14	9.6	8.62	88.4
K3620-03	CL-01-071219-B	2	1.13	9.78	8.98	90.8
K3620-05	CL-01-071219-C	3	1.14	9.62	8.7	89.2
K3620-07	CL-01-071219-D	4	1.13	9.57	8.71	89.8
K3620-09	CL-01-071219-E	5	1.13	9.96	8.97	88.8
K3733-08	KY030LC068-SD-190709	6	1.18	9.98	8.16	79.3
K3733-09	K3733-08MS	7	1.18	9.98	8.16	79.3
K3733-10	K3733-08MSD	8	1.18	9.98	8.16	79.3
K3733-11	KY030LC051-SD-190709	9	1.13	9.73	5.95	56.0
K3733-12	KY030LC001-SD-190709	10	1.17	9.93	7.3	70.0
K3733-13	DUP-0537-190709	11	1.17	9.93	7.3	70.0
K3776-01	FK-SF-01-006-A	12	1.14	9.89	8.85	88.1
K3776-02	FK-SF-01-006-A	13	1.13	9.91	8.9	88.5
K3776-03	FK-SF-01-006-B	14	1.15	9.62	8.52	87.0
K3776-04	FK-SF-01-006-B	15	1.13	9.65	8.66	88.4
K3778-01	STOCKPILE-COMP	16	1.14	9.98	8.85	87.2
K3783-01	KY030LC011-SD-190711	17	1.14	9.94	7.79	75.6
K3783-02	KY030LC015-SD-190711	18	1.14	9.94	7.05	67.2
K3783-03	KY030LC018-SD-190711	19	1.15	9.58	5.43	50.8
K3784-01	OILY-DEBRIS	20	1.00	2.00	2.00	100.0
K3792-01	K580-B5	21	1.17	9.75	8.97	90.9
K3792-02	K580-B6	22	1.17	9.77	8.93	90.2
K3797-05	SVOC-GPC-BLANK	23	1.00	2.00	2.00	100.0
K3797-06	PEST-GPC-BLANK	24	1.00	2.00	2.00	100.0
K3797-07	PEST-GPC-BLANK-SPIKE	25	1.00	2.00	2.00	100.0
K3797-08	PCB-GPC-BLANK	26	1.00	2.00	2.00	100.0
K3797-09	PCB-GPC-BLANK-SPIKE	27	1.00	2.00	2.00	100.0
K3797-10	SVOC-GPC2-BLANK	28	1.00	2.00	2.00	100.0
K3797-11	PEST-GPC2-BLANK	29	1.00	2.00	2.00	100.0
K3797-12	PEST-GPC2-BLANK-SPIKE	30	1.00	2.00	2.00	100.0
K3797-13	PCB-GPC2-BLANK	31	1.00	2.00	2.00	100.0
K3797-14	PCB-GPC2-BLANK -SPIKE	32	1.00	2.00	2.00	100.0
K3801-01	S12-11-B	33	1.15	9.62	8.75	89.7
K3801-02	S12-11-S	34	1.11	9.95	8.96	88.8
K3801-03	S12-12-B	35	1.15	9.56	8.69	89.7
K3801-04	S12-12-S	36	1.13	9.62	8.53	87.2
K3801-05	S12-13-B	37	1.13	9.69	8.79	89.5
K3801-06	S12-13-S	38	1.13	9.98	9.2	91.2
K3801-07	S12-14-B	39	1.13	9.88	9.06	90.6



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:25
In Date: 07/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 07/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB103945

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3801-08	S12-14-S	40	1.13	9.9	9.06	90.4
K3801-09	S12-15-B	41	1.13	9.68	9.02	92.3
K3801-10	S12-15-S	42	1.12	9.58	8.63	88.8
K3801-11	S12-16-B	43	1.14	9.52	8.78	91.2
K3801-12	S12-16-S	44	1.13	9.75	8.82	89.2
K3801-13	S12-17-S	45	1.14	9.89	8.89	88.6
K3801-14	S12-18-S	46	1.16	9.5	8.64	89.7
K3801-15	S12-19-S	47	1.14	9.54	8.87	92.0
K3801-16	S12-20-S	48	1.14	9.85	9.22	92.8
K3801-17	S12-21-S	49	1.13	9.54	9.02	93.8
K3801-18	S12-22-S	50	1.15	9.64	9.00	92.5
K3801-19	S13-01-B	51	1.15	9.97	9.03	89.3
K3801-20	S13-02-B	52	1.14	9.97	9.03	89.4
K3801-21	S13-03-B	53	1.13	9.7	8.84	90.0
K3801-22	S13-04-B	54	1.14	9.6	8.51	87.1
K3801-23	DUP-14	55	1.14	9.63	8.78	90.0
K3801-24	DUP-16	56	1.13	9.91	8.93	88.8
K3802-01	S12-08-S	57	1.13	9.98	8.58	84.2
K3802-02	S12-08-B	58	1.13	9.54	8.51	87.8
K3802-03	S12-03-S-01	59	1.14	9.96	9.33	92.9
K3807-01	FK-PS-01-053	60	1.16	9.96	8.25	80.6
K3807-02	FK-PS-01-053	61	1.13	9.79	7.88	77.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$