



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 07/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 07/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB103962

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3760-06	DB8Y6	1	1.15	9.96	8.04	78.2
K3760-07	DB8Y7	2	1.15	9.68	9.05	92.6
K3760-08	DB8Y8	3	1.2	9.82	8.6	85.8
K3760-09	DB8Y9	4	1.14	9.7	9.17	93.8
K3760-10	DB8Z0	5	1.14	9.72	9.14	93.2
K3760-11	DB8Z1	6	1.16	9.66	6.99	68.6
K3760-12	DB8Z2	7	1.18	9.97	8.39	82.0
K3760-13	DB8Z3	8	1.14	9.87	9.33	93.8
K3760-14	DB8Z4	9	1.14	9.67	8.66	88.2
K3760-15	DB8Z4MS	10	1.14	9.67	8.66	88.2
K3760-16	DB8Z4MSD	11	1.14	9.67	8.66	88.2
K3760-17	DB8Z5	12	1.19	9.98	9.31	92.4
K3760-18	DB8Z6	13	1.18	9.92	9.91	99.9
K3760-19	VHBLK01	14	1.00	2.00	2.00	100.0
K3760-21	DB8W6	15	1.00	2.00	2.00	100.0
K3760-22	DB8W7	16	1.00	2.00	2.00	100.0
K3760-23	DB8W8	17	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$