



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:35
In Date: 07/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 07/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB103964

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3764-01	A52G6	1	1.13	9.89	7.42	71.8
K3764-02	A52G9	2	1.13	9.72	3.19	24.0
K3764-03	A52H5	3	1.15	9.94	5.22	46.3
K3764-04	A52H6	4	1.16	9.97	4.26	35.2
K3764-05	A52H7	5	1.15	9.95	7.64	73.8
K3764-06	A52H8	6	1.13	9.83	7.74	76.0
K3764-07	A52J1	7	1.14	9.9	6.77	64.3
K3764-10	A52S7	8	1.00	2.00	2.00	100.0
K3764-13	A52T1	9	1.00	2.00	2.00	100.0
K3764-15	A52F9	10	1.12	9.82	7.05	68.2
K3764-16	A52G2	11	1.13	9.89	8.29	81.7
K3764-17	A52G4	12	1.18	9.62	7.42	73.9
K3764-18	A52G4MS	13	1.18	9.62	7.42	73.9
K3764-19	A52G4MSD	14	1.18	9.62	7.42	73.9
K3764-20	A52G8	15	1.11	9.92	3.48	26.9
K3764-21	A52H0	16	1.16	9.84	5.03	44.6
K3764-22	A52H1	17	1.18	9.93	7.61	73.5
K3764-23	A52H3	18	1.14	9.98	6.95	65.7
K3764-24	A52H4	19	1.12	9.8	6.33	60.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$