



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/18/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:18
In Date: 07/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 07/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104021

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3615-01	HR-06-071719-A	43	1.11	9.98	9.41	93.6
K3615-03	HR-06-071719-B	44	1.15	9.98	9.48	94.3
K3615-05	HR-06-071719-C	45	1.12	9.94	9.59	96.0
K3618-05	HD-01-071719	46	1.17	9.94	8.74	86.3
K3835-01	SB-1 (0-2)	1	1.15	9.5	8.78	91.4
K3835-02	SB-1 (7-9)	2	1.16	9.96	8.86	87.5
K3835-03	SB-2 (0-2)	3	1.17	9.92	9.02	89.7
K3835-04	SB-2 (10-12)	4	1.13	9.85	8.59	85.6
K3835-05	SB-3 (0-2)	5	1.12	9.88	8.69	86.4
K3835-06	SB-3 (6-8)	6	1.15	9.52	8.38	86.4
K3835-07	SB-4 (0-2)	7	1.14	9.68	8.74	89.0
K3835-08	SB-4 (7-9)	8	1.15	9.94	8.93	88.5
K3835-09	SB-5 (0-2)	9	1.17	9.94	8.84	87.5
K3835-10	SB-5 (7-9)	10	1.16	9.97	8.71	85.7
K3835-23	SB-12 (0-2)	11	1.15	9.88	8.63	85.7
K3835-24	SB-12 (6-8)	12	1.12	9.97	8.96	88.6
K3835-25	SB-9 (0-2)	13	1.15	9.96	9.29	92.4
K3835-26	SB-9 (5-7)	14	1.11	9.7	8.47	85.7
K3835-27	SB-10 (0-2)	15	1.13	9.94	9.25	92.2
K3835-28	SB-10 (8-10)	16	1.13	9.73	8.71	88.1
K3835-29	SB-11 (0-2)	17	1.14	9.97	9.11	90.3
K3835-30	SB-11 (8-10)	18	1.13	9.82	8.35	83.1
K3835-31	SB-7 (0-2)	19	1.1	9.82	9.24	93.3
K3835-32	SB-7 (13-15)	20	1.14	9.55	8.4	86.3
K3861-01	S10-01-S	21	1.12	9.77	9.11	92.4
K3861-02	S10-01-B	22	1.12	9.99	9.3	92.2
K3861-03	S10-02-S	23	1.13	9.58	8.58	88.2
K3861-04	S10-02-B	24	1.13	9.95	8.16	79.7
K3861-05	S15-01-S	25	1.12	9.69	8.38	84.7
K3861-06	S15-01-B	26	1.13	9.9	8.95	89.2
K3862-01	S3-06-B	27	1.13	9.92	9.58	96.1
K3862-02	S3-07-B	28	1.13	9.94	9.39	93.8
K3862-03	S3-08-B	29	1.14	9.63	8.76	89.8
K3862-04	S3-08-S	30	1.13	9.5	9.11	95.3
K3862-05	S3-09-B	31	1.13	9.67	9.04	92.6
K3862-06	S3-09-S	32	1.14	9.83	9.52	96.4
K3862-07	S3-10-B	33	1.12	9.79	9.32	94.6
K3862-08	S3-10-S	34	1.13	9.65	9.16	94.2
K3862-09	S3-11-S	35	1.12	9.95	9.61	96.1



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/18/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:18
In Date: 07/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 07/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104021

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3862-10	S3-12-S	36	1.13	9.5	8.88	92.6
K3862-11	DUP-2	37	1.13	9.95	9.46	94.4
K3862-12	DUP-3	38	1.13	9.83	9.48	96.0
K3866-01	RT-1488	47	1.13	9.69	8.67	88.1
K3867-01	390711744	39	1.12	9.98	2.96	20.8
K3868-01	ASPHALT-1	48	1.00	2.00	2.00	100.0
K3868-02	ASPHALT-2	49	1.00	2.00	2.00	100.0
K3868-03	ASPHALT-3	50	1.00	2.00	2.00	100.0
K3868-04	ASPHALT-4	51	1.00	2.00	2.00	100.0
K3868-05	TP-1	52	1.13	9.51	9.16	95.8
K3868-06	TP-1	53	1.13	9.51	9.16	95.8
K3868-07	TP-2	54	1.14	9.76	9.57	97.8
K3868-08	TP-2	55	1.14	9.76	9.57	97.8
K3868-09	TP-3	56	1.13	9.98	9.08	89.8
K3868-10	TP-3	57	1.13	9.98	9.08	89.8
K3868-11	TP-4	58	1.13	9.87	9.33	93.8
K3868-12	TP-4	59	1.13	9.87	9.33	93.8
K3869-01	390712098	40	1.12	9.9	1.42	3.4
K3871-01	390710709.2	42	1.13	9.97	1.35	2.5
K3872-01	390711720	41	1.13	9.97	1.76	7.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$