



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/19/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 07/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 07/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104045

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3851-01	MBENP0	1	1.12	9.8	9.02	91.0
K3851-02	MBENP1	2	1.16	9.62	7.08	70.0
K3851-03	MBENP1D	3	1.16	9.62	7.08	70.0
K3851-04	MBENP1S	4	1.16	9.62	7.08	70.0
K3851-05	MBENP4	5	1.15	9.59	8.33	85.1
K3851-06	MBENP5	6	1.16	9.6	8.52	87.2
K3851-07	MBENQ0	7	1.17	9.58	6.99	69.2
K3851-08	MBENQ1	8	1.18	9.89	7.09	67.9
K3851-09	MBENQ3	9	1.16	9.96	7.26	69.3
K3851-10	MBENR6	10	1.17	9.8	4.39	37.3
K3851-11	MBENR7	11	1.18	9.97	4.86	41.9
K3851-12	MBENR9	12	1.1	9.92	6.24	58.3
K3851-13	MBENS0	13	1.15	9.7	7.87	78.6
K3851-14	MBENS1	14	1.15	9.92	7.93	77.3
K3851-15	MBENS2	15	1.12	9.59	7.02	69.7
K3851-16	MBENR8	16	1.13	9.63	7.1	70.2
K3851-17	MBENS3	17	1.13	9.68	8.55	86.8
K3851-18	MBENS4	18	1.16	9.98	8.31	81.1
K3851-19	MBENS5	19	1.16	9.7	7.21	70.8
K3851-20	MBENS6	20	1.11	9.54	7.93	80.9
K3851-21	MBENS7	21	1.18	9.93	8.39	82.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$