



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/19/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 07/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 07/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104055

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3850-01	BENP0	1	1.12	9.8	9.02	91.0
K3850-02	BENP1	2	1.16	9.62	7.08	70.0
K3850-03	BENP1MS	3	1.16	9.62	7.08	70.0
K3850-04	BENP1MSD	4	1.16	9.62	7.08	70.0
K3850-05	BENP4	5	1.15	9.59	8.33	85.1
K3850-06	BENP5	6	1.16	9.6	8.52	87.2
K3850-07	BENQ0	7	1.17	9.58	6.99	69.2
K3850-08	BENQ1	8	1.18	9.89	7.09	67.9
K3850-09	BENQ3	9	1.16	9.96	7.26	69.3
K3850-10	BENR6	10	1.17	9.8	4.39	37.3
K3850-11	BENR7	11	1.18	9.97	4.86	41.9
K3850-12	BENR9	12	1.1	9.92	6.24	58.3
K3850-13	BENS0	13	1.15	9.7	7.87	78.6
K3850-14	BENS1	14	1.15	9.92	7.93	77.3
K3850-15	BENS2	15	1.12	9.59	7.02	69.7
K3850-16	BENR8	16	1.13	9.63	7.1	70.2
K3850-17	BENS3	17	1.13	9.68	8.55	86.8
K3850-18	BENS4	18	1.16	9.98	8.31	81.1
K3850-19	BENS5	19	1.16	9.7	7.21	70.8
K3850-20	BENS6	20	1.11	9.54	7.93	80.9
K3850-21	BENS7	21	1.18	9.93	8.39	82.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$