



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/7/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 07/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 07:40
Out Date: 07/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104074

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3864-01	MA52F6	1	1.14	9.71	6.43	61.7
K3864-02	MA52F7	2	1.15	9.82	7.68	75.3
K3864-03	MA52F8	3	1.15	9.89	7.91	77.3
K3864-04	MA52G0	4	1.14	9.66	7.88	79.1
K3864-05	MA52G1	5	1.15	9.96	6.71	63.1
K3864-06	MA52G3	6	1.12	9.59	7.67	77.3
K3864-07	MA52G7	7	1.16	9.71	6.74	65.3
K3864-08	MA52H2	8	1.14	9.84	4.87	42.9
K3864-09	MA52J0	9	1.16	9.97	4.63	39.4
K3864-10	MA52F2	10	1.18	9.97	7.14	67.8
K3864-11	MA52F3	11	1.13	9.93	2.4	14.4
K3864-12	MA52F4	12	1.13	9.96	7.77	75.2
K3864-13	MA52F5	13	1.14	9.98	7.12	67.6
K3864-14	MA52G5	14	1.16	9.56	4.6	41.0
K3864-15	MA52H9	15	1.17	9.82	8.06	79.7
K3864-16	MA52H9D	16	1.17	9.82	8.06	79.7
K3864-17	MA52H9S	17	1.17	9.82	8.06	79.7
K3864-18	MA52T2	18	1.00	2.00	2.00	100.0
K3864-19	MA52T4	19	1.00	2.00	2.00	100.0
K3864-20	MA52T7	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$