



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:05
In Date: 07/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 1.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:00
Out Date: 07/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104118

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3909-01	MH11P8	1	1.00	2.00	2.00	100.0
K3909-02	MH11P9	2	1.00	2.00	2.00	100.0
K3909-03	MH11Q0	3	1.00	2.00	2.00	100.0
K3909-04	MH11Q1	4	1.00	2.00	2.00	100.0
K3909-05	MH11Q2	5	1.00	2.00	2.00	100.0
K3909-06	MH11Q3	6	1.00	2.00	2.00	100.0
K3909-07	MH11Q4	7	1.00	2.00	2.00	100.0
K3909-08	MH11Q5	8	1.00	2.00	2.00	100.0
K3909-09	MH11Q6	9	1.00	2.00	2.00	100.0
K3909-10	MH11Q7	10	1.00	2.00	2.00	100.0
K3909-11	MH11Q8	11	1.00	2.00	2.00	100.0
K3909-12	MH11Q9	12	1.00	2.00	2.00	100.0
K3909-13	MH11R0	13	1.00	2.00	2.00	100.0
K3909-14	MH11R1	14	1.00	2.00	2.00	100.0
K3909-15	MH11R2	15	1.00	2.00	2.00	100.0
K3909-16	MH11R3	16	1.00	2.00	2.00	100.0
K3909-17	MH11R4	17	1.00	2.00	2.00	100.0
K3909-18	MH11R5	18	1.00	2.00	2.00	100.0
K3909-19	MH11R6	19	1.00	2.00	2.00	100.0
K3909-20	MH11R7	20	1.00	2.00	2.00	100.0
K3909-21	MH11R7D	21	1.00	2.00	2.00	100.0
K3909-22	MH11R7S	22	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$