



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:15
In Date: 07/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:05
Out Date: 07/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104119

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3910-01	MH11R8	1	1.00	2.00	2.00	100.0
K3910-02	MH11R9	2	1.00	2.00	2.00	100.0
K3910-03	MH11S0	3	1.00	2.00	2.00	100.0
K3910-04	MH11S1	4	1.00	2.00	2.00	100.0
K3910-05	MH11S2	5	1.00	2.00	2.00	100.0
K3910-06	MH11S3	6	1.00	2.00	2.00	100.0
K3910-07	MH11S4	7	1.00	2.00	2.00	100.0
K3910-08	MH11S5	8	1.00	2.00	2.00	100.0
K3910-09	MH11S6	9	1.00	2.00	2.00	100.0
K3910-10	MH11S7	10	1.00	2.00	2.00	100.0
K3910-11	MH11S8	11	1.00	2.00	2.00	100.0
K3910-12	MH11S9	12	1.00	2.00	2.00	100.0
K3910-13	MH11T0	13	1.00	2.00	2.00	100.0
K3910-14	MH11T1	14	1.00	2.00	2.00	100.0
K3910-15	MH11T2	15	1.00	2.00	2.00	100.0
K3910-16	MH11T2D	16	1.00	2.00	2.00	100.0
K3910-17	MH11T2S	17	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$