



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/23/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 07/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 09:30
Out Date: 07/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104137

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3895-01	MBENY7	1	1.11	9.96	5.81	53.1
K3895-02	MBEP00	2	1.13	9.73	5.45	50.2
K3895-03	MBEP01	3	1.13	9.87	7.03	67.5
K3895-04	MBEP02	4	1.19	9.97	7.81	75.4
K3895-05	MBEP11	5	1.13	9.77	6.15	58.1
K3895-06	MBEP11D	6	1.13	9.77	6.15	58.1
K3895-07	MBEP11S	7	1.13	9.77	6.15	58.1
K3895-08	MBEP12	8	1.14	9.8	5.77	53.5
K3895-09	MBEP13	9	1.12	9.57	5.33	49.8
K3895-10	MBEP14	10	1.13	9.72	6.24	59.5
K3895-11	MBEP15	11	1.1	9.57	6.11	59.1
K3895-12	MBEP16	12	1.13	9.8	5.87	54.7
K3895-13	MBEP17	13	1.17	9.98	6.38	59.1
K3895-14	MBEP18	14	1.15	9.77	6.56	62.8
K3895-15	MBEP19	15	1.11	9.66	8.41	85.4
K3895-16	MBEP20	16	1.12	9.89	8.68	86.2
K3895-17	MBEP21	17	1.11	9.67	6.03	57.5
K3895-18	MBEP22	18	1.17	9.94	3.83	30.3
K3895-19	MBEP23	19	1.18	9.96	3.89	30.9
K3895-20	MBEP24	20	1.17	9.98	4.98	43.2
K3895-21	MBEP32	21	1.13	9.56	5.03	46.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$