



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 7/25/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 15:10

In Date: 07/24/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105

Time OUT: 08:15

Out Date: 07/24/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB104188

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3981-01	ETNE1	1	1.12	9.81	8.93	89.9
K3981-02	ETNE1MS	2	1.12	9.81	8.93	89.9
K3981-03	ETNE1MSD	3	1.12	9.81	8.93	89.9
K3981-04	ETNE2	4	1.13	9.71	9.07	92.5
K3981-05	ETNE3	5	1.14	9.98	9.47	94.2
K3981-06	ETNE4	6	1.15	9.84	9.18	92.4
K3981-07	ETNE5	7	1.16	9.96	9.33	92.8
K3981-08	ETNE6	8	1.14	9.87	8.64	85.9
K3981-09	ETNE7	9	1.13	9.78	8.83	89.0
K3981-10	ETNE8	10	1.13	9.86	8.96	89.7
K3981-11	ETNE9	11	1.14	9.98	9.09	89.9
K3981-12	ETNF0	12	1.15	9.89	8.8	87.5
K3981-13	ETNF1	13	1.14	9.85	9.06	90.9
K3981-14	ETNF2	14	1.16	9.74	8.88	90.0
K3981-15	ETNF3	15	1.13	9.94	9.06	90.0
K3981-16	ETNF4	16	1.13	9.91	8.8	87.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$