



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/26/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 07/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 07/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104220

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4009-01	MH11Z0	1	1.00	2.00	2.00	100.0
K4009-02	MH11Z1	2	1.00	2.00	2.00	100.0
K4009-03	MH11Z2	3	1.00	2.00	2.00	100.0
K4009-04	MH11Z2D	4	1.00	2.00	2.00	100.0
K4009-05	MH11Z2S	5	1.00	2.00	2.00	100.0
K4009-06	MH11Z3	6	1.00	2.00	2.00	100.0
K4009-07	MH11Z4	7	1.00	2.00	2.00	100.0
K4009-08	MH11Z5	8	1.00	2.00	2.00	100.0
K4009-09	MH11Z6	9	1.00	2.00	2.00	100.0
K4009-10	MH11Z7	10	1.00	2.00	2.00	100.0
K4009-11	MH11Z8	11	1.00	2.00	2.00	100.0
K4009-12	MH11Z9	12	1.00	2.00	2.00	100.0
K4009-13	MH1200	13	1.00	2.00	2.00	100.0
K4009-14	MH1201	14	1.00	2.00	2.00	100.0
K4009-15	MH1202	15	1.00	2.00	2.00	100.0
K4009-16	MH1203	16	1.00	2.00	2.00	100.0
K4009-17	MH1204	17	1.00	2.00	2.00	100.0
K4009-18	MH1205	18	1.00	2.00	2.00	100.0
K4009-19	MH1206	19	1.00	2.00	2.00	100.0
K4009-20	MH1207	20	1.00	2.00	2.00	100.0
K4009-21	MH1208	21	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$