



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:38
In Date: 07/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 07/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104245

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3979-01	DB908	1	1.00	2.00	2.00	100.0
K3979-02	DB922	2	1.15	9.94	8.99	89.2
K3979-03	DB923	3	1.12	9.6	8.76	90.1
K3979-04	DB924	4	1.18	9.56	8.52	87.6
K3979-05	DB925	5	1.14	9.84	8.92	89.4
K3979-06	DB926	6	1.14	9.7	8.89	90.5
K3979-07	DB927	7	1.13	9.64	8.53	87.0
K3979-08	DB928	8	1.14	9.73	8.33	83.7
K3979-09	DB929	9	1.12	9.61	8.15	82.8
K3979-10	DB930	10	1.14	9.98	8.96	88.5
K3979-11	DB931	11	1.14	9.73	8.67	87.7
K3979-12	DB932	12	1.13	9.94	8.61	84.9
K3979-13	DB933	13	1.13	9.6	8.67	89.0
K3979-14	DB934	14	1.12	9.54	8.66	89.5
K3979-15	DB935	15	1.13	9.67	8.5	86.3
K3979-16	VHBLK01	16	1.00	2.00	2.00	100.0
K3979-17	DB936	17	1.15	9.85	9.54	96.4
K3979-18	DB937	18	1.13	9.87	9.16	91.9
K3979-19	DB938	19	1.14	9.64	8.4	85.4
K3979-20	DB939	20	1.12	9.63	8.86	91.0
K3979-21	DB940	21	1.13	9.76	9.09	92.2
K3979-22	DB940MS	22	1.13	9.76	9.09	92.2
K3979-23	DB940MSD	23	1.13	9.76	9.09	92.2
K3979-24	DB941	24	1.14	9.56	8.64	89.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$