



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 07/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 07/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104259

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4013-01	TANK-ABC-SLUDGE	1	1.18	9.97	1.89	8.1
K4015-05	SVOC-GPC-BLANK	2	1.00	2.00	2.00	100.0
K4015-06	PEST-GPC-BLANK	3	1.00	2.00	2.00	100.0
K4015-07	PEST-GPC-BLANK-SPIKE	4	1.00	2.00	2.00	100.0
K4015-08	PCB-GPC-BLANK	5	1.00	2.00	2.00	100.0
K4015-09	PCB-GPC-BLANK-SPIKE	6	1.00	2.00	2.00	100.0
K4015-10	SVOC-GPC2-BLANK	7	1.00	2.00	2.00	100.0
K4015-11	PEST-GPC2-BLANK	8	1.00	2.00	2.00	100.0
K4015-12	PEST-GPC2-BLANK-SPIKE	9	1.00	2.00	2.00	100.0
K4015-13	PCB-GPC2-BLANK	10	1.00	2.00	2.00	100.0
K4015-14	PCB-GPC2-BLANK -SPIKE	11	1.00	2.00	2.00	100.0
K4020-01	71701	12	1.00	2.00	2.00	100.0
K4021-01	TP-10	13	1.12	9.71	8.6	87.1
K4021-02	TP-10 (1-1.5)	14	1.15	9.96	8.96	88.6
K4021-03	TP-10 (3.5-4)	15	1.16	9.96	9.11	90.3
K4021-05	TP-9	16	1.14	9.69	8.42	85.1
K4021-06	TP-9 (1-1.5)	17	1.18	9.88	8.97	89.5
K4021-07	TP-9 (3.5-4)	18	1.12	9.6	8.56	87.7
K4021-09	TP-7	19	1.15	9.54	8.4	86.4
K4021-10	TP-7 (1-1.5)	20	1.18	9.83	8.84	88.6
K4021-11	TP-7 (3.5-4)	21	1.1	9.76	8.27	82.8
K4021-13	TP-8-A	22	1.16	9.79	8.74	87.8
K4021-14	TP-8-A (1-1.5)	23	1.11	9.72	8.54	86.3
K4021-15	TP-8-A (3.5-4)	24	1.19	9.81	8.78	88.1
K4024-01	TP-8	25	1.13	9.85	8.84	88.4
K4024-02	TP-8 (1-1.5)	26	1.14	9.94	8.93	88.5
K4024-03	TP-8 (3.5-4)	27	1.13	9.6	8.54	87.5
K4024-05	TP-6-A	28	1.13	9.78	8.67	87.2
K4024-06	TP-6-A (1-1.5)	29	1.15	9.97	9.31	92.5
K4024-07	TP-6-A (3.5-4)	30	1.18	9.6	8.48	86.7
K4027-01	S01-01-B	31	1.13	9.71	8.12	81.5
K4027-02	S01-01-S	32	1.16	9.71	8.12	81.4
K4027-03	DUP-1	33	1.14	9.54	7.79	79.2
K4032-01	S6-05-B	34	1.19	9.69	8.83	89.9
K4032-02	S14-05-S	35	1.15	9.53	8.73	90.5
K4032-03	S14-06-B	36	1.15	9.97	9.25	91.8
K4032-04	S14-06-S	37	1.13	9.73	9.06	92.2
K4032-05	S14-07-B	38	1.15	9.66	8.92	91.3
K4032-06	S14-08-B	39	1.1	9.82	8.82	88.5



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 7/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 07/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 07/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104259

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4032-07	S14-07-S	40	1.15	9.61	8.74	89.7
K4032-08	S14-09-B	41	1.18	9.75	8.92	90.3
K4032-09	S8-03-S	42	1.16	9.96	8.86	87.5
K4032-10	S8-03-B	43	1.14	9.75	8.72	88.0
K4032-11	DUP-8	44	1.12	9.62	8.86	91.1
K4032-12	S12-01-B-01	45	1.15	9.94	9.17	91.2
K4032-13	S9-04-S-01	46	1.17	9.58	8.62	88.6
K4032-14	S4-03-B-01	47	1.15	9.97	9.02	89.2
K4032-15	DUP-11	48	1.18	9.88	8.77	87.2
K4034-01	FK-PS-01-054	49	1.15	9.69	7.95	79.6
K4034-02	FK-PS-01-054	50	1.14	9.81	7.89	77.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$