



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 07/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:33
Out Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104277

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4030-01	MBENZ9	1	1.14	9.74	6.87	66.6
K4030-02	MBEP68	2	1.15	9.76	8.06	80.3
K4030-03	MBEP69	3	1.18	9.79	8.42	84.1
K4030-04	MBENQ6	4	1.18	9.74	7.71	76.3
K4030-05	MBENQ7	5	1.12	9.55	8.06	82.3
K4030-06	MBENY8	6	1.15	9.92	6.59	62.0
K4030-07	MBENY9	7	1.14	9.71	6.88	67.0
K4030-08	MBENZ0	8	1.18	9.77	6.25	59.0
K4030-09	MBENZ1	9	1.16	9.52	5.68	54.1
K4030-10	MBENZ2	10	1.19	9.79	6.37	60.2
K4030-11	MBENZ2D	11	1.19	9.79	6.37	60.2
K4030-12	MBENZ2S	12	1.19	9.79	6.37	60.2
K4030-13	MBENZ3	13	1.15	9.77	6.08	57.2
K4030-14	MBENQ9	14	1.15	9.65	7.88	79.2
K4030-15	MBENR0	15	1.16	9.64	6.01	57.2
K4030-16	MBENR1	16	1.17	9.53	7.1	70.9
K4030-17	MBENR4	17	1.16	9.97	8.48	83.1
K4030-18	MBENR5	18	1.15	9.89	7.4	71.5
K4030-19	MBEP71	19	1.12	9.63	6.28	60.6
K4030-20	MBENW3	20	1.14	9.6	5.44	50.8
K4030-21	MBENW4	21	1.11	9.58	5.38	50.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$