



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:10
In Date: 07/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104285

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4031-01	MBENQ8	1	1.14	9.7	7.14	70.1
K4031-02	MBENQ8D	2	1.14	9.7	7.14	70.1
K4031-03	MBENQ8S	3	1.14	9.7	7.14	70.1
K4031-04	MBENX5	4	1.13	9.69	6.04	57.4
K4031-05	MBENX6	5	1.12	9.55	7.09	70.8
K4031-06	MBENX7	6	1.15	9.95	8.17	79.8
K4031-07	MBENX8	7	1.13	9.95	5.84	53.4
K4031-08	MBENX9	8	1.17	9.93	8.26	80.9
K4031-09	MBENY0	9	1.12	9.56	7.79	79.0
K4031-10	MBENX4	10	1.1	9.73	7.86	78.3
K4031-11	MBENY1	11	1.16	9.98	8.41	82.2
K4031-12	MBENY2	12	1.14	9.82	8.43	84.0
K4031-13	MBENY3	13	1.19	9.69	6.85	66.6
K4031-14	MBENY4	14	1.1	9.62	7.58	76.1
K4031-15	MBENY5	15	1.15	9.76	6.06	57.0
K4031-16	MBENY6	16	1.14	9.59	5.31	49.3
K4031-17	MBENT6	17	1.13	9.6	5.95	56.9
K4031-18	MBENT7	18	1.15	9.82	4.63	40.1
K4031-19	MBENT8	19	1.17	9.98	5.88	53.5
K4031-20	MBENT9	20	1.11	9.75	6.06	57.3
K4031-21	MBENW0	21	1.16	9.76	6.29	59.7
K4031-22	MBENW2	22	1.13	9.88	8.75	87.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$