



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 07/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104293

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3621-07	OK-01-072919	28	1.12	9.88	9.44	95.0
K3622-13	OR-03-072919-A	29	1.14	9.74	9.3	94.9
K3622-15	OR-03-072919-B	30	1.13	9.88	9.34	93.8
K4044-01	TP-6 (1-1.5)	1	1.14	9.86	9.06	90.8
K4044-02	TP-6 (3.5-4)	2	1.13	9.91	9.32	93.3
K4044-03	TP-6	3	1.13	9.87	8.84	88.2
K4044-05	TP-5	4	1.13	9.54	8.76	90.7
K4044-07	TP-5 (1-1.5)	5	1.14	9.96	9.48	94.6
K4044-08	TP-5 (3.5-4)	6	1.14	9.7	8.23	82.8
K4044-09	TP-1	7	1.14	9.82	9.22	93.1
K4044-11	TP-1 (1-1.5)	8	1.13	9.8	9.14	92.4
K4044-12	TP-1 (3.5-4)	9	1.14	9.97	9.24	91.7
K4044-13	TP-4 (1-1.5)	10	1.14	9.5	8.48	87.8
K4044-14	TP-4 (3.5-4)	11	1.13	9.61	7.67	77.1
K4044-15	TP-4	12	1.13	9.76	8.57	86.2
K4044-17	TP-3	13	1.13	9.7	8.55	86.6
K4044-19	TP-3 (1-1.5)	14	1.13	9.54	8.85	91.8
K4044-20	TP-3 (3.5-4)	15	1.13	9.98	8.92	88.0
K4046-01	S1-02-B	31	1.12	9.55	8.05	82.2
K4046-02	S1-02-S	32	1.15	9.6	8.03	81.4
K4046-03	S1-03-S	33	1.16	9.94	9.00	89.3
K4046-04	S2-01-B	34	1.14	9.96	8.29	81.1
K4046-05	S2-01-S	35	1.18	9.71	7.98	79.7
K4046-06	DUP-2	36	1.15	9.5	7.86	80.4
K4046-07	S2-02-S	37	1.15	9.74	7.79	77.3
K4046-08	S3-01-B	38	1.13	9.85	8.22	81.3
K4046-09	S3-01-S	39	1.15	9.92	8.4	82.7
K4046-10	S3-02-B	40	1.13	9.81	8.65	86.6
K4046-11	S3-02-S	41	1.1	9.82	8.7	87.2
K4047-01	S4-01-B	42	1.15	9.96	9.37	93.3
K4047-02	S4-02-B	43	1.18	9.97	9.61	95.9
K4047-03	S3-07-B-01	44	1.18	9.96	8.37	81.9
K4049-01	SP-1	16	1.14	9.86	8.25	81.5
K4049-02	SP-2	17	1.13	9.56	8.05	82.1
K4049-03	SP-3	18	1.13	9.95	8.46	83.1
K4049-04	SP-4	19	1.14	9.86	8.48	84.2
K4049-05	SP-5	20	1.13	9.55	8.41	86.5
K4049-06	SP-6	21	1.14	9.76	8.64	87.0
K4049-07	SP-7	22	1.14	9.63	8.29	84.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/30/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:20
In Date: 07/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104293

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4049-08	SP-8	23	1.14	9.85	8.83	88.3
K4049-09	SP-9	24	1.13	9.63	8.44	86.0
K4049-10	SP-10	25	1.14	9.98	8.71	85.6
K4049-11	SP-11	26	1.13	9.75	8.52	85.7
K4049-12	SP-12	27	1.12	9.98	8.66	85.1
K4050-24	COMP-4	45	1.00	2.00	2.00	100.0
K4051-05	COMP-1	46	1.00	2.00	2.00	100.0
K4051-16	COMP-2	47	1.00	2.00	2.00	100.0
K4051-23	COMP-3	48	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$