



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/31/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:48
In Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 07/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104303

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K3618-09	HD-02-073019	60	1.13	9.63	8.88	91.2
K4043-13	BP-VPB203-GW-258-260	1	1.15	9.98	3.58	27.5
K4043-14	BP-VPB203-GW-298-300	2	1.18	9.9	3.42	25.7
K4043-15	BP-VPB203-GW-418-420	3	1.14	9.97	2.67	17.3
K4052-01	S4-01-A	4	1.12	9.66	7.96	80.1
K4052-02	S4-01-B	5	1.16	9.72	7.42	73.1
K4052-03	S4-02-B	6	1.14	9.93	7.51	72.5
K4052-06	DUP-4	7	1.18	9.56	7.69	77.7
K4053-01	S12-01-A	8	1.14	9.59	6.48	63.2
K4053-02	S12-02-A	9	1.15	9.66	7.44	73.9
K4054-01	S6-01-A	10	1.18	9.85	7.7	75.2
K4055-01	S10-01-A	11	1.14	9.89	9.22	92.3
K4055-03	DUP-10	12	1.17	9.5	8.87	92.4
K4056-01	S7-01-A	13	1.15	9.54	7.06	70.4
K4057-01	S13-02-A	14	1.15	9.79	7.88	77.9
K4057-02	S13-01-B	15	1.17	9.81	8.21	81.5
K4057-03	S13-02-B	16	1.11	9.53	7.83	79.8
K4057-04	S13-01-C	17	1.17	9.77	8.76	88.3
K4057-05	S13-02-C	18	1.18	9.51	8.06	82.6
K4057-06	S13-03-C	19	1.14	9.94	8.65	85.3
K4057-07	S13-04-C	20	1.1	9.57	7.74	78.4
K4058-01	S9-01-A	21	1.15	9.58	8.79	90.6
K4058-02	S9-02-A	22	1.13	9.71	9.41	96.5
K4058-03	S9-03-A	23	1.14	9.77	9.1	92.2
K4058-04	S9-04-A	24	1.16	9.57	8.92	92.3
K4058-05	S9-05-A	25	1.14	9.53	8.77	90.9
K4059-01	S5-01-A	26	1.17	9.69	5.22	47.5
K4059-02	S5-02-A	27	1.14	9.52	6.4	62.8
K4059-03	S5-03-A	28	1.18	9.58	5.76	54.5
K4060-01	S8-01-A	29	1.14	9.56	7.99	81.4
K4061-01	S1-01-A	30	1.19	9.6	9.02	93.1
K4061-02	S1-02-A	31	1.15	9.77	8.98	90.8
K4061-03	S1-03-A	32	1.1	9.58	8.97	92.8
K4061-04	S1-04-A	33	1.14	9.59	7.89	79.9
K4061-05	S1-05-A	34	1.11	9.52	8.01	82.0
K4061-16	DUP-1	35	1.18	9.73	9.08	92.4
K4062-01	S2-01-A	36	1.13	9.76	8.14	81.2
K4062-02	S2-02-A	37	1.15	9.7	8.3	83.6
K4062-04	S2-04-A	38	1.14	9.5	7.01	70.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/31/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:48
In Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 07/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104303

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4062-07	S2-01-B	39	1.16	9.87	8.82	87.9
K4062-08	S2-02-B	40	1.13	9.55	8.15	83.4
K4062-10	S2-04-B	41	1.12	9.92	6.58	62.0
K4062-11	S2-05-B	42	1.18	9.75	8.16	81.4
K4062-12	S2-06-B	43	1.15	9.72	7.89	78.6
K4062-13	S2-01-C	44	1.18	9.6	8.35	85.2
K4062-14	S2-02-C	45	1.15	9.73	8.69	87.9
K4062-15	S2-03-C	46	1.18	9.55	7.92	80.5
K4062-16	S2-04-C	47	1.14	9.88	7.83	76.5
K4062-17	S2-05-C	48	1.16	9.64	7.75	77.7
K4062-18	S2-06-C	49	1.14	9.85	7.79	76.3
K4062-19	DUP-2	50	1.12	9.51	7.78	79.4
K4063-01	S3-01-B	51	1.16	9.58	7.5	75.3
K4063-02	S3-02-B	52	1.15	9.96	7.94	77.1
K4064-01	S11-01-A	53	1.16	9.82	8.44	84.1
K4064-02	S11-02-A	54	1.17	9.8	9.26	93.7
K4064-03	S11-01-B	55	1.15	9.5	8.37	86.5
K4064-04	S11-02-B	56	1.17	9.76	8.73	88.0
K4064-05	S11-03-B	57	1.15	9.96	9.36	93.2
K4064-06	S11-04-B	58	1.11	9.66	8.24	83.4
K4064-07	S11-05-B	59	1.19	9.61	7.79	78.4
K4066-01	TP-1-A	61	1.13	9.57	8.79	90.8
K4066-03	TP-1-A(1-1.5)	62	1.13	9.85	9.06	90.9
K4066-04	TP-1-A(3.5-4)	63	1.14	9.81	8.93	89.9
K4066-05	TP-11(1.0-1.5)	64	1.13	9.87	9.02	90.3
K4066-07	TP-11	65	1.14	9.6	8.55	87.6
K4066-10	TP-12(3.5-4.0)	66	1.13	9.98	8.92	88.0
K4066-11	TP-12	67	1.13	9.98	9.2	91.2
K4066-13	TP-13(1.0-1.5)	68	1.12	9.75	8.63	87.0
K4066-14	TP-13(3.5-4.0)	69	1.16	9.67	8.48	86.0
K4066-15	TP-13	70	1.15	9.96	7.52	72.3
K4072-01	S4-01-S	71	1.14	9.68	7.48	74.2
K4072-02	S4-01-B	72	1.14	9.78	8.98	90.7
K4072-03	S4-02-S	73	1.14	9.71	7.53	74.6
K4074-01	SS-01_072919	74	1.15	9.88	9.16	91.8
K4074-02	SS-02_072919	75	1.15	9.97	9.09	90.0
K4074-03	SS-03_072919	76	1.15	9.96	8.93	88.3
K4074-04	SS-04_072919	77	1.14	9.83	8.95	89.9
K4074-05	SS-05_072919	78	1.14	9.74	9.11	92.7



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/31/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:48
In Date: 07/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 07/31/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104303

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4074-06	SS-06_072919	79	1.14	9.94	8.69	85.8
K4074-07	SS-07_072919	80	1.13	9.99	9.04	89.3
K4074-08	SS-08_072919	81	1.13	9.98	9.6	95.7
K4074-13	SS-13_073019	82	1.14	9.6	8.59	88.1
K4074-14	SS-14_073019	83	1.12	9.97	8.98	88.8
K4076-01	S2-01-B	84	1.12	9.98	8.18	79.7
K4076-02	S2-02-B	85	1.13	9.58	8.16	83.2
K4077-02	20190612-LOC-13	86	1.13	9.59	8.78	90.4
K4077-06	20190612-LOC-13	87	1.13	9.59	8.78	90.4
K4078-01	INFO-BOX-4-073019	88	1.00	2.00	2.00	100.0
K4078-02	INFO-BOX-5-073019	89	1.00	2.00	2.00	100.0
K4078-03	INFO-BOX-6-073019	90	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$