



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/2/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:35
In Date: 08/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104370

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4116-01	JLLD0	1	1.13	9.7	7.49	74.2
K4116-02	JLLD1	2	1.14	9.83	7.77	76.3
K4116-03	JLLD2	3	1.16	9.84	7.42	72.1
K4116-04	JLLD3	4	1.14	9.81	7.7	75.7
K4116-05	JLLD4	5	1.14	9.88	7.68	74.8
K4116-06	JLLD5	6	1.12	9.74	7.6	75.2
K4116-07	JLLD6	7	1.14	9.89	7.82	76.3
K4116-08	JLLD7	8	1.13	9.72	7.97	79.6
K4116-09	JLLD9	9	1.13	9.82	7.5	73.3
K4116-10	JLLE0	10	1.13	9.94	7.89	76.7
K4116-11	JLLE1	11	1.13	9.75	7.68	76.0
K4116-12	JLLE4	12	1.14	9.94	6.49	60.8
K4116-13	JLLE5	13	1.14	9.72	7.76	77.2
K4116-14	JLLE6	14	1.13	9.72	7.79	77.5
K4116-15	JLLE7	15	1.13	9.85	7.89	77.5
K4116-16	JLLE8	16	1.12	9.84	8.02	79.1
K4116-17	JLLE9	17	1.14	9.92	7.65	74.1
K4116-18	JLLF0	18	1.14	9.76	7.71	76.2
K4116-21	VHBLK01	19	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$