



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 08:20
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 14:55
Out Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104383

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4137-01	S15-10-B	1	1.18	9.72	9.03	91.9
K4137-02	S15-09-B	2	1.18	9.63	8.76	89.7
K4137-03	S15-07-S	3	1.16	9.5	8.46	87.5
K4137-04	S15-08-B	4	1.18	9.6	8.72	89.5
K4137-05	S15-06-S	5	1.18	9.68	8.74	88.9
K4137-06	S15-07-B	6	1.15	9.54	8.7	90.0
K4137-07	15-06-B	7	1.18	9.64	8.83	90.4
K4137-08	S15-05-S	8	1.17	9.77	9.19	93.3
K4137-09	S15-04-S	9	1.16	9.94	9.32	92.9
K4137-10	S15-05-B	10	1.12	9.68	8.77	89.4
K4137-11	S15-03-S	11	1.17	9.51	8.92	92.9
K4137-12	S15-02-S	12	1.16	9.77	9.25	94.0
K4137-13	S15-04-B	13	1.14	9.98	9.29	92.2
K4137-14	S11-03-B-02	14	1.16	9.96	9.21	91.5
K4137-15	DUP-11	15	1.17	9.93	8.99	89.3
K4137-16	S5-01-B	16	1.14	9.56	9.22	96.0
K4137-17	S5-01-S	17	1.14	9.7	9.42	96.7
K4137-18	S5-02-B	18	1.11	9.9	9.03	90.1
K4137-19	S5-04-B	19	1.15	9.5	8.95	93.4
K4137-20	S5-03-B	20	1.19	9.57	8.5	87.2
K4137-21	S5-02-S	21	1.18	9.54	9.15	95.3
K4138-01	S6-01-B	22	1.16	9.94	8.99	89.2
K4138-02	S6-01-S	23	1.18	9.67	8.00	80.3
K4138-03	S6-02-B	24	1.13	9.95	8.87	87.8
K4138-04	S6-02-S	25	1.14	9.51	7.88	80.5
K4138-05	DUP-6	26	1.18	9.66	8.57	87.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$