



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4090-01	OR-02-080219-A	218	1.13	9.69	9.3	95.4
K4090-03	OR-02-080219-B	219	1.14	9.63	8.94	91.9
K4095-01	SU-02-080219-A	220	1.14	9.87	9.06	90.7
K4095-03	SU-02-080219-B	221	1.13	9.71	8.71	88.3
K4129-01	1S-H- (0-1)	1	1.13	9.82	8.85	88.8
K4129-02	1S-H- (30)	2	1.14	9.65	8.57	87.3
K4129-03	2S-E1-E4- (0-1) COMP	3	1.13	9.79	8.79	88.5
K4129-04	2S-E1- (0-1)	4	1.14	9.81	8.93	89.9
K4129-05	2S-E2- (0-1)	5	1.13	9.65	8.6	87.7
K4129-06	2S-E3- (0-1)	6	1.12	9.75	8.91	90.3
K4129-07	2S-E4- (0-1)	7	1.12	9.91	8.99	89.5
K4129-08	2S-E1-E4- (30) COMP	8	1.13	9.9	8.78	87.2
K4129-09	2S-E1- (30)	9	1.13	9.96	9.18	91.2
K4129-10	2S-E2- (30)	10	1.13	9.74	8.55	86.2
K4129-11	2S-E3- (30)	11	1.11	9.98	8.94	88.3
K4129-12	2S-E4- (30)	12	1.11	9.88	8.98	89.7
K4130-01	1S-F1-F4-COMP- (0-1)	13	1.12	9.77	8.64	86.9
K4130-02	1S-F1- (0-1)	14	1.11	9.86	8.9	89.0
K4130-03	1S-F2- (0-1)	15	1.12	9.79	8.51	85.2
K4130-04	1S-F3- (0-1)	16	1.11	9.9	9.04	90.2
K4130-05	1S-F4- (0-1)	17	1.12	9.93	8.79	87.1
K4130-06	1S-F1-F4-COMP- (30)	18	1.13	9.94	8.82	87.3
K4130-07	1S-F1- (30)	19	1.13	9.74	8.77	88.7
K4130-08	1S-F2- (30)	20	1.14	9.9	8.69	86.2
K4130-09	1S-F3- (30)	21	1.13	9.9	8.88	88.4
K4130-10	1S-F4- (30)	22	1.14	9.86	8.68	86.5
K4130-11	1S-G- (0-30) -COMP	23	1.12	9.92	8.87	88.1
K4130-12	1S-G- (0-1)	24	1.13	9.72	8.61	87.1
K4130-13	1S-G- (30)	25	1.15	9.83	8.83	88.5
K4130-14	1S-L- (0-30) -COMP	26	1.13	9.9	8.85	88.0
K4130-15	1S-L- (0-1)	27	1.14	9.62	8.51	86.9
K4130-16	1S-L- (30)	28	1.14	9.69	8.59	87.1
K4130-17	1S-M- (0-30) -COMP	29	1.15	9.68	8.46	85.7
K4130-18	1S-M- (0-1)	30	1.13	9.96	8.83	87.2
K4130-19	1S-M- (30)	31	1.14	9.84	8.43	83.8
K4130-20	1S-H- (0-30) -COMP	32	1.13	9.56	8.47	87.1
K4131-01	1S-D1-D4-COMP- (0-1)	33	1.14	9.76	8.81	89.0
K4131-02	1S-D1- (0-1)	34	1.13	9.81	8.93	89.9
K4131-03	1S-D2- (0-1)	35	1.13	9.91	8.82	87.6



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4131-04	1S-D3- (0-1)	36	1.14	9.79	8.9	89.7
K4131-05	1S-D4- (0-1)	37	1.13	9.96	8.8	86.9
K4131-06	1S-D1-D4-COMP- (30)	38	1.13	9.78	8.81	88.8
K4131-07	1S-D1- (30)	39	1.13	9.75	8.87	89.8
K4131-08	1S-D2- (30)	40	1.15	9.97	9.01	89.1
K4131-09	1S-D3- (30)	41	1.13	9.98	9.02	89.2
K4131-10	1S-D4- (30)	42	1.14	9.67	8.81	89.9
K4131-11	1S-E1-E4-COMP- (0-1)	43	1.16	9.95	9.11	90.4
K4131-12	1S-E1- (0-1)	44	1.14	9.75	9.00	91.3
K4131-13	1S-E2- (0-1)	45	1.14	9.7	8.79	89.4
K4131-14	1S-E3- (0-1)	46	1.13	9.98	9.17	90.8
K4131-15	1S-E4- (0-1)	47	1.14	9.55	8.4	86.3
K4131-16	1S-E1-E4-COMP- (30)	48	1.13	9.55	8.68	89.7
K4131-17	1S-E1- (30)	49	1.12	9.77	8.92	90.2
K4131-18	1S-E2- (30)	50	1.15	9.84	8.91	89.3
K4131-19	1S-E3- (30)	51	1.13	9.69	8.72	88.7
K4131-20	1S-E4- (30)	52	1.14	9.59	8.72	89.7
K4132-01	BOX-6-080119-5	53	1.00	2.00	2.00	100.0
K4132-02	BOX-6-080119-9	54	1.00	2.00	2.00	100.0
K4132-03	BOX-6-080119-11	55	1.00	2.00	2.00	100.0
K4132-04	BOX-6-080119-34	56	1.00	2.00	2.00	100.0
K4132-05	BOX-6-080119-39	57	1.00	2.00	2.00	100.0
K4132-06	BOX-6-080119-52	58	1.00	2.00	2.00	100.0
K4133-01	CLAM-BUCKET-080119	59	1.00	2.00	2.00	100.0
K4133-02	MINI-EX-080119	60	1.00	2.00	2.00	100.0
K4133-03	DEWATERING-SCREEN-1-080119	61	1.00	2.00	2.00	100.0
K4134-01	BOX-5-080119-5	62	1.00	2.00	2.00	100.0
K4134-02	BOX-5-080119-11	63	1.00	2.00	2.00	100.0
K4134-03	BOX-5-080119-19	64	1.00	2.00	2.00	100.0
K4134-04	BOX-5-080119-36	65	1.00	2.00	2.00	100.0
K4134-05	BOX-5-080119-45	66	1.00	2.00	2.00	100.0
K4134-06	BOX-5-080119-53	67	1.00	2.00	2.00	100.0
K4135-01	1S-A1-A4-COMP- (0-1)	68	1.13	9.5	8.56	88.8
K4135-02	1S-A1- (0-1)	69	1.13	9.64	8.63	88.1
K4135-03	1S-A2- (0-1)	70	1.12	9.93	9.01	89.6
K4135-04	1S-A3- (0-1)	71	1.13	9.51	8.27	85.2
K4135-05	1S-A4- (0-1)	72	1.13	9.55	8.7	89.9
K4135-06	1S-A1-A4-COMP- (30)	73	1.13	9.67	8.77	89.5
K4135-07	1S-A1- (30)	74	1.15	9.98	9.28	92.1

**PERCENT SOLID**

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4135-08	1S-A2- (30)	75	1.13	9.7	8.68	88.1
K4135-09	1S-A3- (30)	76	1.13	9.5	8.34	86.1
K4135-10	1S-A4- (30)	77	1.13	9.57	8.74	90.2
K4135-11	1S-B1-B4-COMP- (0-1)	78	1.14	9.51	8.71	90.4
K4135-12	1S-B1- (0-1)	79	1.13	9.61	8.51	87.0
K4135-13	1S-B2- (0-1)	80	1.13	9.73	8.92	90.6
K4135-14	1S-B3- (0-1)	81	1.11	9.55	8.42	86.6
K4135-15	1S-B4- (0-1)	82	1.14	9.88	8.53	84.6
K4135-16	1S-B1-B4-COMP- (30)	83	1.13	9.59	8.79	90.5
K4135-17	1S-B1- (30)	84	1.14	9.78	8.71	87.6
K4135-18	1S-B2- (30)	85	1.14	9.81	9.09	91.7
K4135-19	1S-B3- (30)	86	1.14	9.59	8.56	87.8
K4135-20	1S-B4- (30)	87	1.14	9.92	8.74	86.6
K4136-01	1S-C1-C4-COMP- (0-1)	88	1.13	9.65	8.88	91.0
K4136-02	1S-C1- (0-1)	89	1.12	9.85	8.96	89.8
K4136-03	1S-C2- (0-1)	90	1.14	9.98	9.11	90.2
K4136-04	1S-C3- (0-1)	91	1.12	9.96	8.93	88.3
K4136-05	1S-C4- (0-1)	92	1.14	9.73	9.32	95.2
K4136-06	1S-C1-C4-COMP- (30)	93	1.12	9.77	8.82	89.0
K4136-07	1S-C1- (30)	94	1.12	9.51	8.66	89.9
K4136-08	1S-C2- (30)	95	1.12	9.53	8.42	86.8
K4136-09	1S-C3- (30)	96	1.14	9.73	8.73	88.4
K4136-10	1S-C4- (30)	97	1.15	9.81	9.21	93.1
K4139-05	SVOC-GPC-BLANK	131	1.00	2.00	2.00	100.0
K4139-06	PEST-GPC-BLANK	132	1.00	2.00	2.00	100.0
K4139-07	PEST-GPC-BLANK-SPIKE	133	1.00	2.00	2.00	100.0
K4139-08	PCB-GPC-BLANK	134	1.00	2.00	2.00	100.0
K4139-09	PCB-GPC-BLANK-SPIKE	135	1.00	2.00	2.00	100.0
K4139-10	SVOC-GPC2-BLANK	136	1.00	2.00	2.00	100.0
K4139-11	PEST-GPC2-BLANK	137	1.00	2.00	2.00	100.0
K4139-12	PEST-GPC2-BLANK-SPIKE	138	1.00	2.00	2.00	100.0
K4139-13	PCB-GPC2-BLANK	139	1.00	2.00	2.00	100.0
K4139-14	PCB-GPC2-BLANK -SPIKE	140	1.00	2.00	2.00	100.0
K4140-21	P001-SS008-1824-02	176	1.14	9.81	9.41	95.4
K4140-22	P001-SS015-0006-01	177	1.13	9.97	8.55	83.9
K4140-23	P001-SS015-2430-01	178	1.13	9.65	9.16	94.2
K4140-24	P001-TP2-1824-02	179	1.14	9.65	5.58	52.2
K4140-25	P001-TP2-4248-01	180	1.14	9.66	8.02	80.8
K4140-26	P001-TP3-0006-01	181	1.13	9.73	2.98	21.5



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4140-27	P001-TP3-0624-01	182	1.13	9.55	4.33	38.0
K4140-28	P001-TP5-1224-01	183	1.14	9.66	6.83	66.8
K4140-29	P001-TP5-2436-01	184	1.13	9.62	3.39	26.6
K4140-30	P001-TP6-0108-01	185	1.13	9.85	8.6	85.7
K4140-31	P001-TP6-2026-01	186	1.13	9.98	5.5	49.4
K4140-32	P001-TP8-0006-01	187	1.14	9.62	9.24	95.5
K4140-33	P001-TP8-1824-01	188	1.13	9.97	9.48	94.5
K4140-34	P001-SS001-0006-01	189	1.14	9.74	9.22	94.0
K4140-35	P001-SS001-1824-01	190	1.14	9.75	9.47	96.7
K4140-36	P001-SS002-0006-01	191	1.13	9.7	9.38	96.3
K4140-37	P001-SS002-1824-01	192	1.14	9.96	9.67	96.7
K4140-38	P001-SS003-0006-01	193	1.13	9.98	9.79	97.9
K4140-39	P001-SS003-1824-01	194	1.13	9.96	9.55	95.4
K4140-40	P001-SS04-0006-01	195	1.13	9.69	8.95	91.4
K4141-01	P001-SS004-1824-01	196	1.12	9.95	9.18	91.3
K4141-02	P001-SS005-0006-01	197	1.13	9.91	7.22	69.4
K4141-03	P001-SS005-1824-01	198	1.12	9.86	8.88	88.8
K4141-04	P001-SS006-0006-01	199	1.13	9.8	7.99	79.1
K4141-05	P001-SS006-1824-01	200	1.13	9.98	7.49	71.9
K4141-06	P001-SS007-0006-01	201	1.14	9.78	9.36	95.1
K4141-07	P001-SS007-1824-01	202	1.14	9.69	9.27	95.1
K4141-08	P001-SS008-0006-01	203	1.15	9.72	9.16	93.5
K4141-09	P001-SS008-1824-01	204	1.14	9.54	8.85	91.8
K4141-10	P001-SS009-0006-01	205	1.14	9.92	9.65	96.9
K4141-11	P001-SS009-1218-01	206	1.16	9.81	9.11	91.9
K4141-12	P001-SS010-0006-01	207	1.16	9.68	9.19	94.2
K4141-13	P001-SS010-1824-01	208	1.15	9.65	9.15	94.1
K4141-14	P001-SS011-0006-01	209	1.15	9.78	9.04	91.4
K4141-15	P001-SS011-1824-01	210	1.16	9.81	9.51	96.5
K4141-16	P001-SS012-0006-01	211	1.15	9.73	8.96	91.0
K4141-17	P001-SS012-1824-01	212	1.14	9.74	7.49	73.8
K4141-18	P001-SS013-0006-01	213	1.15	9.69	7.88	78.8
K4141-19	P001-SS013-1824-01	214	1.15	9.82	9.32	94.2
K4141-20	P001-SS014-0006-01	215	1.13	9.81	9.14	92.3
K4141-21	P001-SS014-1824-01	216	1.12	9.87	9.45	95.2
K4141-22	P001-TP2-1824-01	217	1.14	9.56	5.22	48.5
K4143-01	PST-006-SW125-073019	98	1.12	9.72	9.28	94.9
K4143-02	PST-006-SW125-073019-1	99	1.16	9.9	9.36	93.8
K4143-03	PST-006-SW126-073019	100	1.13	9.55	9.01	93.6



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4143-04	PST-006-SW127-073019	101	1.14	9.63	9.12	94.0
K4143-05	PST-006-B259-073019	102	1.18	9.8	9.44	95.8
K4143-06	PST-006-B260-073019	103	1.14	9.71	9.28	95.0
K4143-07	PST-006-B261-073019	104	1.13	9.5	9.00	94.0
K4143-08	EXT-006-SW130-073019	105	1.13	9.54	8.77	90.8
K4143-09	EXT-006-SW131-073019	106	1.13	9.67	9.01	92.3
K4143-10	EXT-006-SW132-073019	107	1.13	9.64	8.83	90.5
K4143-11	PST-006-SW052-073019	108	1.15	9.89	9.42	94.6
K4143-12	PST-006-SW053-073019	109	1.14	9.67	9.18	94.3
K4143-13	EXT-003-SW075-073019	110	1.15	9.83	9.26	93.4
K4143-14	EXT-003-SW074-073019	111	1.15	9.71	9.22	94.3
K4143-15	EXT-004-SW073-073019	112	1.14	9.62	9.16	94.6
K4143-16	EXT-004-SW072-073019	113	1.15	9.73	9.13	93.0
K4143-17	EXT-023-SW071-073019	114	1.15	9.65	9.22	94.9
K4143-18	EXT-023-SW076-073019	115	1.14	9.59	9.3	96.6
K4143-19	EXT-008-SW118-073119	116	1.11	9.75	9.22	93.9
K4143-20	EXT-008-SW119-073119	117	1.18	9.96	9.48	94.5
K4144-01	EXT-006-SW129-073119	118	1.19	9.61	9.08	93.7
K4144-02	EXT-006-SW127-073119	119	1.11	9.72	9.44	96.7
K4144-03	EXT-007-SW125-073119	120	1.15	9.94	9.55	95.6
K4144-04	PST-006-B258-073119	121	1.16	9.97	9.42	93.8
K4144-05	EXT-006-SW128-073119	122	1.16	9.7	9.26	94.8
K4144-06	EXT-008-SW121-073119	123	1.17	9.97	9.58	95.6
K4144-07	EXT-008-SW120-073119	124	1.13	9.75	9.16	93.2
K4144-08	EXT-008-SW115-073119	125	1.13	9.64	9.32	96.2
K4144-09	PST-008-B250-073119	126	1.17	9.58	8.67	89.2
K4144-10	PST-008-SW120-073119	127	1.14	9.8	9.12	92.1
K4144-11	EXT-008-SW111-073119	128	1.15	9.74	8.74	88.4
K4144-12	EXT-008-SW109-073119	129	1.14	9.97	9.3	92.4
K4144-13	EXT-008-SW108-073119	130	1.14	9.65	9.22	94.9
K4145-08	BP-VPB203-GW-618-620	141	1.14	9.98	2.57	16.2
K4145-09	BP-VPB203-GW-638-640	142	1.14	9.97	2.7	17.7
K4147-01	C0AM4	143	1.13	9.67	7.29	72.1
K4147-02	C0AM5	144	1.13	9.56	7.53	75.9
K4147-03	C0AM5MS	145	1.13	9.56	7.53	75.9
K4147-04	C0AM5MSD	146	1.13	9.56	7.53	75.9
K4147-05	C0AM6	147	1.13	9.85	7.09	68.3
K4147-06	C0AM8	148	1.13	9.58	6.77	66.7
K4147-07	C0AM9	149	1.13	9.68	5.9	55.8



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4147-08	C0AN0	150	1.14	9.61	6.36	61.6
K4147-09	C0AN1	151	1.13	9.55	7.7	78.0
K4150-01	S15-01-B	152	1.14	9.96	9.39	93.5
K4150-02	S15-02-B	153	1.12	9.5	8.52	88.3
K4150-03	DUP-18	154	1.13	9.91	9.36	93.7
K4150-04	S15-03-B	155	1.15	9.66	8.79	89.8
K4150-05	S15-01-S	156	1.18	9.68	9.21	94.5
K4150-06	S17-01-B	157	1.15	9.7	8.8	89.5
K4150-07	S17-02-B	158	1.12	9.7	8.82	89.7
K4150-08	S17-01-S	159	1.11	9.66	8.44	85.7
K4150-09	S17-04-B	160	1.11	9.71	9.26	94.8
K4150-10	DUP-4	161	1.11	9.92	9.26	92.5
K4150-11	S17-03-B	162	1.18	9.79	8.96	90.4
K4150-12	S17-02-S	163	1.1	9.63	9.09	93.7
K4150-13	S8-01-S-01	164	1.15	9.5	9.15	95.8
K4151-01	S5-01-S	165	1.15	9.51	7.8	79.5
K4151-02	S5-01-B	166	1.14	9.81	8.68	87.0
K4151-03	S5-02-S	167	1.15	9.71	7.8	77.7
K4151-04	S5-02-B	168	1.15	9.74	8.72	88.1
K4151-05	DUP-4	169	1.14	9.56	8.6	88.6
K4151-06	S5-03-S	170	1.13	9.81	8.52	85.1
K4151-07	S5-04-B	171	1.14	9.82	8.44	84.1
K4151-08	S5-04-S	172	1.15	9.67	7.83	78.4
K4151-09	DUP-5	173	1.14	9.53	7.47	75.4
K4151-10	S5-05-S	174	1.14	9.58	8.35	85.4
K4151-11	S5-05-B	175	1.16	9.61	8.06	81.7
K4153-01	RT-4221	222	1.14	9.72	8.22	82.5
K4153-02	RT-4221	223	1.14	9.72	8.22	82.5
K4155-01	TS-01_080219	224	1.00	2.00	2.00	100.0
K4156-01	RT1494	225	1.14	9.98	8.15	79.3
K4156-02	RT1494	226	1.14	9.98	8.15	79.3
K4156-04	RT1873	227	1.13	9.78	8.45	84.6
K4156-05	RT1873	228	1.13	9.78	8.45	84.6
K4156-07	RT3380	229	1.15	9.84	7.93	78.0
K4156-08	RT3380	230	1.15	9.84	7.93	78.0
K4156-10	RT3456	231	1.13	9.89	7.64	74.3
K4156-11	RT3456	232	1.13	9.89	7.64	74.3
K4156-13	RBR-200035	233	1.14	9.73	8.98	91.3
K4156-14	RBR-200035	234	1.14	9.73	8.98	91.3



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/5/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 08/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:40
Out Date: 08/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: % SOLIDS-OVEN

QC:LB104389

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K4158-01	FK-PS-01-055	235	1.13	9.65	7.96	80.2
K4158-02	FK-PS-01-055	236	1.14	9.98	8.31	81.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$