



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/6/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 08/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104418

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4042-01	MBENW1	1	1.15	9.84	7.95	78.3
K4042-02	MBENW1D	2	1.15	9.84	7.95	78.3
K4042-03	MBENW1S	3	1.15	9.84	7.95	78.3
K4042-04	MBENS8	4	1.16	9.86	3.87	31.1
K4042-05	MBENS9	5	1.17	9.97	4.09	33.2
K4042-06	MBENT0	6	1.12	9.83	4.32	36.7
K4042-07	MBENT1	7	1.18	9.73	3.88	31.6
K4042-08	MBENT2	8	1.14	9.83	3.93	32.1
K4042-09	MBENT3	9	1.15	9.77	4.26	36.1
K4042-10	MBENT4	10	1.14	9.51	4.77	43.4
K4042-11	MBENT5	11	1.15	9.74	4.52	39.2
K4042-12	MBENW5	12	1.14	9.53	8.42	86.8
K4042-13	MBENW6	13	1.18	9.65	6.88	67.3
K4042-14	MBENW7	14	1.15	9.5	7.66	78.0
K4042-15	MBENW8	15	1.14	9.54	7.89	80.4
K4042-16	MBENW9	16	1.15	9.96	7.97	77.4
K4042-17	MBENX0	17	1.16	9.97	8.3	81.0
K4042-18	MBENX1	18	1.12	9.66	8.54	86.9
K4042-19	MBENX2	19	1.16	9.95	5.25	46.5
K4042-20	MBENX3	20	1.18	9.77	4.8	42.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$