

**PERCENT SOLID**

Supervisor: apatel
Analyst: jignesh
Date: 8/6/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:09
In Date: 08/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104421

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4152-01	2S-H1-H4-COMP- (0-1)	1	1.14	9.94	8.9	88.2
K4152-02	2S-H1- (0-1)	2	1.12	9.66	8.58	87.4
K4152-03	2S-H2- (0-1)	3	1.13	9.77	8.5	85.3
K4152-04	2S-H3- (0-1)	4	1.18	9.79	8.64	86.6
K4152-05	2S-H4- (0-1)	5	1.17	9.7	8.38	84.5
K4152-06	2S-H1-H4-COMP- (30)	6	1.11	9.83	8.71	87.2
K4152-07	2S-H1- (30)	7	1.18	9.9	8.86	88.1
K4152-08	2S-H2- (30)	8	1.16	9.78	8.54	85.6
K4152-09	2S-H3- (30)	9	1.17	9.65	8.57	87.3
K4152-10	2S-H4- (30)	10	1.14	9.68	8.53	86.5
K4152-11	2S-E1-E4-COMP- (0-1)	11	1.16	9.78	8.88	89.6
K4152-12	2S-E1- (0-1)	12	1.15	9.94	9.05	89.9
K4152-13	2S-E2- (0-1)	13	1.13	9.57	8.99	93.1
K4152-14	2S-E3- (0-1)	14	1.17	9.57	8.86	91.5
K4152-15	2S-E4- (0-1)	15	1.13	9.81	9.03	91.0
K4152-16	2S-E1-E4-COMP- (30)	16	1.14	9.74	9.12	92.8
K4152-17	2S-E1- (30)	17	1.1	9.77	9.29	94.5
K4152-18	2S-E2- (30)	18	1.14	9.5	8.69	90.3
K4152-19	2S-E3- (30)	19	1.12	9.78	8.9	89.8
K4152-20	2S-E4- (30)	20	1.11	9.68	9.03	92.4
K4154-01	PS-1	21	1.18	9.89	4.16	34.2
K4154-02	PS-2	22	1.15	9.98	5.83	53.0
K4154-03	PS-3	23	1.16	9.94	6.24	57.9
K4154-04	PS-4	24	1.14	9.96	8.47	83.1
K4154-05	PS-5	25	1.13	9.55	5.55	52.5
K4154-06	PS-6	26	1.17	9.63	7.31	72.6
K4154-07	PS-7	27	1.14	9.75	6.78	65.5
K4154-08	PS-8	28	1.15	9.98	8.34	81.4
K4154-09	PS-9	29	1.16	9.95	6.14	56.7
K4154-10	PS-10	30	1.17	9.6	5.36	49.7
K4162-01	RT-3407	31	1.13	9.71	8.88	90.3
K4162-02	RT-3407	32	1.13	9.71	8.88	90.3
K4162-04	RBR-2000050	33	1.12	9.59	7.69	77.6
K4162-05	RBR-2000050	34	1.12	9.59	7.69	77.6
K4164-01	RBR-2000030	35	1.17	9.78	8.72	87.7
K4164-02	RBR-2000030	36	1.17	9.78	8.72	87.7
K4164-04	RBR-2000051	37	1.19	9.98	9.24	91.6
K4164-05	RBR-2000051	38	1.19	9.98	9.24	91.6
K4169-01	MJ-GB02	39	1.14	9.55	8.5	87.5



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/6/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:09
In Date: 08/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:27
Out Date: 08/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104421

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4169-02	MJ-GB01	40	1.19	9.97	7.89	76.3
K4170-01	S7-01-B	41	1.15	9.83	9.32	94.1
K4170-02	S7-02-B	42	1.14	9.59	9.16	94.9
K4170-03	S7-03-B	43	1.13	9.66	9.26	95.3
K4170-04	S7-04-B	44	1.14	9.57	9.09	94.3
K4170-05	S7-05-B	45	1.15	9.94	9.27	92.4
K4170-06	S7-06-B	46	1.13	9.66	9.39	96.8
K4170-07	S7-07-B	47	1.14	9.76	9.31	94.8
K4170-08	S7-01-S	48	1.14	9.51	9.4	98.7
K4170-09	S7-02-S	49	1.14	9.52	9.44	99.0
K4170-10	S7-03-S	50	1.13	9.67	9.44	97.3
K4170-11	DUP-7	51	1.13	9.86	9.58	96.8
K4170-12	S17-03-S	52	1.14	9.67	9.05	92.7
K4170-13	S17-05-B	53	1.14	9.97	9.37	93.2
K4171-01	S7-01-S	54	1.14	9.82	8.88	89.2
K4171-02	DUP-7	55	1.15	9.58	8.64	88.8
K4171-03	S7-01-B	56	1.15	9.98	9.06	89.6
K4171-04	S7-02-S	57	1.15	9.79	8.89	89.6
K4171-05	S7-02-B	58	1.14	9.67	8.77	89.4
K4171-06	S7-03-S	59	1.14	9.94	8.83	87.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$