



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 14:35
In Date: 08/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:20
Out Date: 08/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104436

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4175-01	MH1216	1	1.00	2.00	2.00	100.0
K4175-02	MH1217	2	1.00	2.00	2.00	100.0
K4175-03	MH1218	3	1.00	2.00	2.00	100.0
K4175-04	MH1219	4	1.00	2.00	2.00	100.0
K4175-05	MH1220	5	1.00	2.00	2.00	100.0
K4175-06	MH1221	6	1.00	2.00	2.00	100.0
K4175-07	MH1222	7	1.00	2.00	2.00	100.0
K4175-08	MH1223	8	1.00	2.00	2.00	100.0
K4175-09	MH1224	9	1.00	2.00	2.00	100.0
K4175-10	MH1225	10	1.00	2.00	2.00	100.0
K4175-11	MH1225D	11	1.00	2.00	2.00	100.0
K4175-12	MH1225S	12	1.00	2.00	2.00	100.0
K4175-13	MH1226	13	1.00	2.00	2.00	100.0
K4175-14	MH1227	14	1.00	2.00	2.00	100.0
K4175-15	MH1228	15	1.00	2.00	2.00	100.0
K4175-16	MH1229	16	1.00	2.00	2.00	100.0
K4175-17	MH1230	17	1.00	2.00	2.00	100.0
K4175-18	MH1231	18	1.00	2.00	2.00	100.0
K4175-19	MH1232	19	1.00	2.00	2.00	100.0
K4175-20	MH1233	20	1.00	2.00	2.00	100.0
K4175-21	MH1234	21	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$