



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/7/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:25
Out Date: 08/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104437

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4176-01	MH1235	1	1.00	2.00	2.00	100.0
K4176-02	MH1236	2	1.00	2.00	2.00	100.0
K4176-03	MH1237	3	1.00	2.00	2.00	100.0
K4176-04	MH1238	4	1.00	2.00	2.00	100.0
K4176-05	MH1239	5	1.00	2.00	2.00	100.0
K4176-06	MH1239D	6	1.00	2.00	2.00	100.0
K4176-07	MH1239S	7	1.00	2.00	2.00	100.0
K4176-08	MH1240	8	1.00	2.00	2.00	100.0
K4176-09	MH1241	9	1.00	2.00	2.00	100.0
K4176-10	MH1242	10	1.00	2.00	2.00	100.0
K4176-11	MH1243	11	1.00	2.00	2.00	100.0
K4176-12	MH1244	12	1.00	2.00	2.00	100.0
K4176-13	MH1245	13	1.00	2.00	2.00	100.0
K4176-14	MH1246	14	1.00	2.00	2.00	100.0
K4176-15	MH1248	15	1.00	2.00	2.00	100.0
K4176-16	MH1249	16	1.00	2.00	2.00	100.0
K4176-17	MH1250	17	1.00	2.00	2.00	100.0
K4176-18	MH1251	18	1.00	2.00	2.00	100.0
K4176-19	MH1252	19	1.00	2.00	2.00	100.0
K4176-20	MH1253	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$