



PERCENT SOLID

Supervisor: apatel  
Analyst: jignesh  
Date: 8/8/2019

OVENTEMP IN Celsius(°C): 108  
Time IN: 16:03  
In Date: 08/07/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:10  
Out Date: 08/08/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %SOLIDS-OVEN

QC:LB104472

| Lab ID   | Client Sample ID | Dish# | Dish Wt(g)<br>(A) | Dish +<br>Sample Wt(g)<br>(B) | Dish + Dry<br>Sample Wt(g)<br>(C) | % Solid |
|----------|------------------|-------|-------------------|-------------------------------|-----------------------------------|---------|
| K4213-01 | ETP63            | 1     | 1.15              | 9.76                          | 8.69                              | 87.6    |
| K4213-02 | ETP64            | 2     | 1.13              | 9.77                          | 8.57                              | 86.1    |
| K4213-03 | ETP65            | 3     | 1.13              | 9.98                          | 8.62                              | 84.6    |
| K4213-04 | ETP66            | 4     | 1.13              | 9.85                          | 8.74                              | 87.3    |
| K4213-05 | ETP66MS          | 5     | 1.13              | 9.85                          | 8.74                              | 87.3    |
| K4213-06 | ETP66MSD         | 6     | 1.13              | 9.85                          | 8.74                              | 87.3    |
| K4213-07 | ETP67            | 7     | 1.12              | 9.53                          | 8.79                              | 91.2    |
| K4213-08 | ETP68            | 8     | 1.14              | 9.67                          | 8.89                              | 90.9    |
| K4213-09 | ETP69            | 9     | 1.13              | 9.87                          | 9.01                              | 90.2    |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$