



PERCENT SOLID

Supervisor: apatel
Analyst: Dodley
Date: 8/12/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:51
In Date: 08/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 09:45
Out Date: 08/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104498

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4183-01	C0AF1	1	1.13	9.87	8.49	84.2
K4183-02	C0AF3	2	1.14	9.58	8.22	83.9
K4183-03	C0AH0	3	1.12	9.64	4.99	45.4
K4183-04	C0AH3	4	1.14	9.98	4.12	33.7
K4183-05	C0AH4	5	1.14	9.72	4.05	33.9
K4183-06	C0AH6	6	1.14	9.54	6.71	66.3
K4183-07	C0AH7	7	1.14	9.89	4.17	34.6
K4183-08	C0AH9	8	1.14	9.98	4.15	34.0
K4183-09	C0AH5	9	1.12	9.63	8.56	87.4
K4183-10	C0AE8	10	1.13	9.69	8.64	87.7
K4183-11	C0AE9	11	1.14	9.81	8.64	86.5
K4183-12	C0AF0	12	1.14	9.69	8.33	84.1
K4183-13	C0AF2	13	1.14	9.71	8.34	84.0
K4183-14	C0AF4	14	1.14	9.82	8.4	83.6
K4183-15	C0AG1	15	1.16	9.9	8.99	89.6
K4183-16	C0AG3	16	1.16	9.8	8.85	89.0
K4183-17	C0AG4	17	1.16	9.81	8.37	83.4
K4183-18	C0AG8	18	1.16	9.65	6.79	66.3
K4183-19	C0AK4	19	1.13	9.9	7.98	78.1
K4183-20	C0AG5	20	1.13	9.73	8.79	89.1
K4183-21	C0AG5MS	21	1.13	9.73	8.79	89.1
K4183-22	C0AG5MSD	22	1.13	9.73	8.79	89.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$