



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 8/28/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:25
In Date: 08/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 08/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: MSC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104500

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
K4216-01	METOB4	1	1.14	9.84	8.48	84.4	
K4216-02	METOB4D	2	1.14	9.84	8.48	84.4	
K4216-03	METOB4S	3	1.14	9.84	8.48	84.4	
K4216-04	METOB5	4	1.12	9.77	7.07	68.8	
K4216-05	METOB6	5	1.12	9.74	9.21	93.9	
K4216-06	METOB7	6	1.12	9.88	8.6	85.4	
K4216-07	METOB8	7	1.12	9.74	8.51	85.7	
K4216-08	METOB9	8	1.12	9.88	8.37	82.8	
K4216-09	METOC0	9	1.12	9.71	8.49	85.8	
K4216-10	METOC1	10	1.15	9.84	8.65	86.3	
K4216-11	METOC2	11	1.13	9.79	8.35	83.4	
K4216-12	METOC3	12	1.12	9.91	9.38	94.0	
K4216-13	METOC4	13	1.1	9.91	8.42	83.1	
K4216-14	METOC8	14	1.11	9.55	8.11	82.9	
K4216-15	METOC9	15	1.12	9.59	7.86	79.6	
K4216-16	METOD0	16	1.12	9.83	7.53	73.6	
K4216-17	METOD2	17	1.00	2.00	2.00	100.0	
K4216-18	METOD3	18	1.11	9.94	4.17	34.7	
K4216-19	METOD5	19	1.12	9.96	4.1	33.7	
K4216-20	METOD6	20	1.11	9.89	7.88	77.1	
K4216-21	METOC5	21	1.13	9.83	8.78	87.9	
K4216-22	METOC6	22	1.13	9.67	8.27	83.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$