



PERCENT SOLID

Supervisor: apatel
Analyst: Dodley
Date: 8/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:33
In Date: 08/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 09:35
Out Date: 08/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104501

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4252-01	ETP80	1	1.14	9.89	8.8	87.5
K4252-02	ETP81	2	1.15	9.73	8.53	86.0
K4252-03	ETP82	3	1.14	9.88	8.53	84.6
K4252-04	ETP83	4	1.15	9.71	8.29	83.4
K4252-05	ETP84	5	1.14	9.71	8.36	84.2
K4252-06	ETP85	6	1.16	9.53	8.27	84.9
K4252-07	ETP85MS	7	1.16	9.53	8.27	84.9
K4252-08	ETP85MSD	8	1.16	9.53	8.27	84.9
K4252-09	ETP86	9	1.15	9.67	8.09	81.5
K4252-10	ETP87	10	1.14	9.91	8.54	84.4
K4252-11	ETP88	11	1.14	9.7	8.66	87.9
K4252-12	ETP89	12	1.13	9.72	8.68	87.9
K4252-13	ETP90	13	1.13	9.79	8.00	79.3
K4252-14	ETP91	14	1.14	9.84	8.9	89.2
K4252-15	ETP92	15	1.14	9.56	7.78	78.9
K4252-16	ETP93	16	1.14	9.58	7.81	79.0
K4252-17	ETP94	17	1.13	9.87	7.96	78.1
K4252-18	ETP95	18	1.14	9.92	8.11	79.4
K4252-19	ETP96	19	1.13	9.75	7.65	75.6
K4252-20	ETP97	20	1.15	9.81	8.72	87.4
K4252-21	ETP98	21	1.14	9.9	8.2	80.6
K4252-22	ETP99	22	1.14	9.84	8.69	86.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$