



PERCENT SOLID

Supervisor: apatel
Analyst: Dodley
Date: 8/12/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:38
In Date: 08/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 10:00
Out Date: 08/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104523

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4286-01	ETPA1	1	1.15	9.72	8.31	83.5
K4286-02	ETPA2	2	1.14	9.91	8.49	83.8
K4286-03	ETPA3	3	1.13	9.6	8.25	84.1
K4286-04	ETPA4	4	1.13	9.79	8.43	84.3
K4286-05	ETPA5	5	1.15	9.62	8.1	82.1
K4286-06	ETPA6	6	1.14	9.92	8.73	86.4
K4286-07	ETPA7	7	1.15	9.9	8.86	88.1
K4286-08	ETPA8	8	1.14	9.69	8.37	84.6
K4286-09	ETPA9	9	1.13	9.83	8.16	80.8
K4286-10	ETPB0	10	1.14	9.65	8.3	84.1
K4286-11	ETPB1	11	1.13	9.78	8.01	79.5
K4286-12	ETPB2	12	1.13	9.75	8.31	83.3
K4286-13	ETPB3	13	1.13	9.87	8.37	82.8
K4286-14	ETPB4	14	1.13	9.79	8.38	83.7
K4286-15	ETPB5	15	1.15	9.98	8.44	82.6
K4286-16	ETPB5MS	16	1.15	9.98	8.44	82.6
K4286-17	ETPB5MSD	17	1.14	9.98	8.44	82.6
K4286-18	ETPB6	18	1.13	9.7	8.8	89.5
K4286-19	ETPB7	19	1.13	9.69	8.72	88.7
K4286-20	ETPB8	20	1.13	9.53	8.74	90.6
K4286-21	ETPB9	21	1.13	9.51	8.78	91.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$