



PERCENT SOLID

Supervisor: apatel
Analyst: Dodley
Date: 8/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 21:32
In Date: 08/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 10:45
Out Date: 08/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB104554

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K4215-01	ETOB4	1	1.14	9.64	8.32	84.5
K4215-02	ETOB4MS	2	1.14	9.64	8.32	84.5
K4215-03	ETOB4MSD	3	1.14	9.64	8.32	84.5
K4215-04	ETOB5	4	1.14	9.51	6.94	69.3
K4215-05	ETOB6	5	1.15	9.95	8.96	88.8
K4215-06	ETOB7	6	1.14	9.62	8.35	85.0
K4215-07	ETOB8	7	1.13	9.87	8.52	84.6
K4215-08	ETOB9	8	1.16	9.76	8.09	80.6
K4215-09	ETOC0	9	1.16	9.98	8.66	85.0
K4215-10	ETOC1	10	1.14	9.51	8.35	86.1
K4215-11	ETOC2	11	1.14	9.73	8.33	83.7
K4215-12	ETOC3	12	1.15	9.61	8.98	92.6
K4215-13	ETOC4	13	1.14	9.56	8.44	86.7
K4215-14	VHBLK01	14	1.00	2.00	2.00	100.0
K4215-15	ETOC8	15	1.15	9.91	8.3	81.6
K4215-16	ETOC9	16	1.15	9.66	7.97	80.1
K4215-17	ETOD0	17	1.14	9.75	7.25	71.0
K4215-18	ETOD3	18	1.14	9.65	4.56	40.2
K4215-19	ETOD5	19	1.15	9.81	4.57	39.5
K4215-20	ETOD6	20	1.15	9.62	8.00	80.9
K4215-21	ETOC5	21	1.14	9.59	8.43	86.3
K4215-22	ETOC6	22	1.15	9.76	8.6	86.5
K4215-23	ETOC7	23	1.15	9.87	8.73	86.9
K4215-24	ETOD4	24	1.16	9.58	8.5	87.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$